

MACINTOSH FIBEROPTIC

- Stainless steel construction with Satin non-glare finish
- Fiberoptic light transmission provides cool, obstruction free viewing
- 2.5v Halogen illumination for true tissue color and long lasting lamp performance
- Removable light pipe facilitates cleaning and service
- Color coded green to signify compatibility ONLY with fiberoptic handles



ORDERING INFORMATION				
Item#	Blade	Description	Dim A	Dim B
4070F	0	Premie	80	55
4071F	1	Infant	92	67
4072F	2	Child	108	83
4073F	3	Med. Adult	126	101
4074F	4	Lrg. Adult	155	134

(Dimensions in millimeters)


MILLER FIBEROPTIC

- Stainless steel construction with Satin non-glare finish
- Fiberoptic light transmission provides cool, obstruction free viewing
- 2.5v Halogen illumination for true tissue color and long lasting lamp performance
- Removable light pipe facilitates cleaning and service
- Color coded green to signify compatibility ONLY with fiberoptic handles



ORDERING INFORMATION				
Item#	Blade	Description	Dim A	Dim B
4080F	0	Premie	78	53
4081F	1	Infant	102	77
4082F	2	Child	155	103
4083F	3	Med. Adult	195	170
4084F	4	Lrg. Adult	205	180

(Dimensions in millimeters)



FIBEROPTIC BATTERY HANDLES

- Light source in handle simplifies blade cleaning, allows for use of brighter, whiter light
- Satin chrome-plated brass handle with a knurled finish for a positive grip
- Color coded green to signify compatibility ONLY with fiberoptic blades
- Available in traditional "C", "AA", and stubby ("AA")
- Batteries not included



ORDERING INFORMATION	
Item#	Description
4065F	Medium "C" size Fiberoptic battery handle
4066F	Small "AA" size Fiberoptic battery handle
4067F	Stubby "AA" size Fiberoptic battery handle
4532-1	Fiberoptic Replacement Lamp, Pkg. of 1

