



ADC COMPLETE 2.5v INSTRUMENT SET

A full featured 2.5v set with 5 different instruments in a fitted case.

The complete set features

- Bayonet locking mounts for each instrument head to allow secure connection to the battery handle
- "C" cell powered, satin finish chrome-plated battery handle with knurled finish
- Rheostatic on/off switch
- Complete in zippered case with fitted foam insert that is scuff and scratch resistant and cleans easily
- Weighs 23½ oz
- Inspected and packaged in the U.S.A. from foreign components
- Batteries NOT included

Otoscope

- Chrome-plated Otoscope head with removable 4x magnifying lens, and insufflator fitting
- Auto-clavable, polypropylene ear specula in 2.5mm, 3.5mm, and 4.5mm sizes
- Adapter to allow use with Welch Allyn or compatible disposable specula
- Nasal speculum adapter with manual screw opening for instrumentation connects to the otoscope head

- Welch Allyn disposable polypropylene ear specula in 2.75mm and 4.25mm sizes are convenient, safe, and economical

Ophthalmoscope

- Ophthalmoscope head with lens wheel featuring 25 corrective lenses from -25 to +40 diopters

Bent Arm Illuminator

- Chrome-plated bent arm larynx and nasopharynx illuminator for lighting the oral cavity and pharyngeal area
- One 3.0mm and one 4.0mm mirror with integral tongue depressor holder that attaches to the illuminator
- One tongue depressor attachment for the illuminator
- One disposable tongue depressor adapter for the illuminator



ORDERING INFORMATION

#5215	Complete 2.5v Diagnostic Instrument Set
-------	---