

moblvac[®] III CS

DESCRIPTION AND APPLICATION

Moblvac III CS is a high vacuum & flow suction unit. It produces constant suction with vacuum & flow surpassing standard wall suction. Moblvac III CS is packaged in a compact, durable molded Polystyrene case. The heart of the Moblvac III CS is a fan cooled carbon rotary vane pump. Moblvac III CS is insulated to provide one of the quietest sounding suction pumps on the market today. It is an excellent choice for the O.R., surgicenters, and doctors offices.

FEATURES AND BENEFITS

- Constant suction providing 550 mmHg and over 38 lpm flow
- Quite operation
- Easy to clean and process for quick turnarounds
- Fan cooled pump, mechanical safety overflow trap and exhaust port filter
- **3 year warranty**

IDEAL FOR

G.I. Labs, Operating Rooms and Surgicenters requiring high vacuum and flow rates with quite performance. Surgical, scoping, and most specialty procedures.



PORTABLE SUCTION SYSTEM

Wall Vacuum Performance

G.I. Labs • Surgi-Centers • Physician's Offices •
O.R., E.R., and/or General Wall Vacuum Backup

TECHNICAL SPECIFICATIONS

Performance

Vacuum Range: 0 to 550mmHg
 Free Air FLOW: 38 LPM minimum
 Pump: fan cooled rotary carbon vane type

Collection Device

Canister: disposable plastic with mechanical shutoff
 Capacity: 1200cc standard
 Tubing: 14" with bacterial filter

Controls

Vacuum Regulator: rotary type on panel
 Vacuum Gauge: Calibrated in mm Hg

Electrical Requirements

AC: 120 volts, 60 Hz, 3A
 220 volts, 50 Hz, 2A

Physical

Weight: 32 lbs
 Dimensions: 28"H x 14"W x 16"D

Available Accessories

- Canisters (disposable & glass)
- Aspirating Change Out Kits
- Filters & tubing
- Dual Vacuum Port

PURCHASING SPECIFICATIONS

Moblvac® IIIcs Part Number: 757000

ORDERING INFORMATION

Order Number:
 6700-1100-900

SHIPPING INFORMATION

Shipping weight:
 50 lbs.

Package Size:
 50" H x 18" W x 18" D

Packaging:
 One per carton

Terms:
 net 30 days

