

AM Series Hi-lo Treatment Tables
share these features:

- 34 oz. heavy duty, hand crafted vinyl upholstery.
- 1-1/2" firm density foam (unless otherwise noted).
- Radiused, seamless corners on upholstered components.
- Footswitch control.
- Power assisted Caster System (unless noted otherwise).
- 120 volts, 60HZ (220 volts, 50 HZ available as special order).
- 400 lb. weight capacity (unless noted otherwise).



The AM-200 Treatment Table

Two Piece Top Section.....	27" x 76"
Minimum Height.....	18"
Maximum Height.....	37"
Width	27"
Length.....	76"
One Piece Head Section.....	15" x 27"
Raises 45 degrees. Lowers 90 degrees	
One Piece Body Section	61" x 27"



AM-100 Treatment Table

One piece top section	27" x 76"
Minimum Height.....	17"
Maximum Height.....	36"
Width	27"
Length.....	76"
Non-skid adjustable glides	



AM-150 Treatment Table

One piece top section	27" x 76"
Minimum Height	18"
Maximum Height	37"
Width	27"
Length	76"