


AM-234 BARIATRIC TABLE

2 Piece Top Section.....34" x 76"
 Minimum Height.....18"
 Maximum Height.....37"
 Width34"
 Length.....76"
 One piece Head / Backrest.....34" x 27".
 Raises 74 degrees.
 Contoured face / nose opening.
 One Piece Body section34" x 49".
 Capacity.....500 lbs.



The AM-250 Treatment Table

Four Piece Top Section.....27" x 76"
 Minimum Height.....18"
 Maximum Height.....37"
 Width27"
 Length.....76"
 Three Piece Head Section.....
 Center Section.....15" x 13"
 Arm Rest.....12 1/2" x 5 1/2"
 Raises 45 degrees. Lowers 90 degrees
 Contoured Face/Nose Opening
 One Piece Body Section61" x 27"



AM-240 BO-BATH TREATMENT TABLE

2 Piece Top Section 40" x 76"
 Minimum Height 18"
 Maximum Height 37"
 Width 40"
 Length 76"
 One piece head/backrest 40"x27"
 Raises 74 degrees
 One Piece Body Section 40"x49"
 Capacity 500 lbs

