

The AM-Series Hi-Lo Treatment Tables

These treatment tables are what put Armedica on the map...strong, stylish and affordable.

For the clinician, the AM-Series offers full height adjustment and section mobility so the patient can be placed in the proper position for maximum results with minimum effort.

For the patient, comfort and easy access getting on and off of the table. All AM-Series Treatment Tables share these features:

- 34 oz. Heavy-duty, handcrafted vinyl upholstery with Permablok³® bacterial protection.
- 1-1/2" Firm density foam tops.
- Radiused, seamless corners on upholstered components.

- Top sections have welded steel support frames.
- Pivoting joints have hardened steel shoulder bolts with locking nuts to ensure years of trouble-free operation.
- Footswitch control.
- Power assisted caster system (unless noted otherwise).
- 110 volts, 60Hz motor (220 volts, 50 Hz available on special order).
- 400 lb. Lifting weight capacity (unless noted otherwise).



The AM-300 Treatment Table

Three Piece Top Section	27" x 76"
Minimum Height	18"
Maximum Height	37"
Width	27"
Length	76"
One Piece Head Section	15" x 27"
Raises 45 degrees. Lowers 90 degrees	
Contoured Face/ Nose Opening.	
One Piece Center Section	22" x 27"
Raises 30 degrees and has a Spring Loaded Lock that must be released before the section will raise.	
One Piece Foot Section	39" x 27"
Raises 70 degrees	