


## AM Series Hi-lo Treatment Tables

### share these features:

- 34 oz. heavy duty, hand crafted vinyl upholstery.
- 1-1/2" firm density foam (unless otherwise noted).
- Radiused, seamless corners on upholstered components.
- Footswitch control.
- Power assisted Caster System (unless noted otherwise).
- 120 volts, 60HZ (220 volts, 50 HZ available as special order).
- 400 lb. weight capacity (unless noted otherwise).

### The AM-350 Treatment Table

Three Piece Top Section.....	27" x 76"
Minimum Height.....	18"
Maximum Height.....	37"
Width.....	27"
Length.....	76"
One Piece Head Section.....	15" x 27"
Raises 45 degrees. Lowers 90 degrees	
Contoured Face/Nose Opening	
One Piece Center Section.....	22" x 27"
Does Not Raise	
One Piece Foot Section.....	39" x 27"
Raises 70 degrees	



### The AM-353 Treatment Table

Three Piece Top Section	27" x 76"
Minimum Height	17"
Maximum Height	36"
Width	27"
Length	76"
One Piece Head Section	15" x 27"
Raises 45 degrees. Lowers 90 degrees	
Contoured Face/Nose Opening	
One Piece Center Section	22" x 27"
Does Not Raise	
One Piece Foot Section	39" x 27"
Raises 70 degrees	
Non-Skid Adjustable Glides	

