

AM Series Hi-lo Treatment Tables

share these features:

- 34 oz. heavy duty, hand crafted vinyl upholstery.
- 1-1/2" firm density foam (unless otherwise noted).
- Radiused, seamless corners on upholstered components.
- Footswitch control.
- Power assisted Caster System (unless noted otherwise).
- 120 volts, 60HZ (220 volts, 50 HZ available as special order).
- 400 lb. weight capacity (unless noted otherwise).

The AM-400 Traction Table

Four Piece Top Section	27 x 76
Minimum Height	18"
Maximum Height	37"
Width	27"
Length Upholstered Top	76"
Length Overall	99"
One Piece Head Section	15" x 27"
Raises 45 degrees. Lowers 90 degrees	
Contoured Face/ Nose Opening.	
Thoracic Section	17-1/2 x 27"
Center Section	10" x 27"
Foot Section	33" x 27"
Raises 80 degrees	



Friction Free Lumbar Section Separates 5" and rolls on 4 Steel Ball Roller Bearings. The Rolling Section can be locked in any position with a positive gear rack system.

The AM-450 Traction Table

Six Piece Top Section	27 x 76
Minimum Height	18"
Maximum Height	37"
Width	27"
Length Upholstered Top	76"
Length Overall	99"
Three Piece Head Section	
Center	15" x 13"
Arm Rest	12-1/2" x 5 1/2"
Raises 45 degrees. Lowers 90 degrees	
Contoured Face/ Nose Opening.	
Thoracic Section	17 1/2 x 27"
Center Section	10" x 27"
Foot Section	33" x 27"
Raises 80 degrees	



Friction Free Lumbar Section Separates 5" and rolls on 4 Steel Ball Roller Bearings. The Rolling Section can be locked in any position with a positive gear rack system.