


## AM SP Series Hi-lo Treatment Tables

include the following features:

- 34 oz. heavy duty vinyl hand crafted upholstery.
- 1-1/2" firm density foam (unless noted otherwise)
- Radiused, seamless corners on upholstered components.
- Footswitch control
- 400 lb. lifting/weight capacity
- 120 volts, 60HZ (220 volts, 50HZ available as special order).
- One-step Caster System




### AM-SP100 TREATMENT TABLE

One piece top section .....	27" x 76"
Minimum Height.....	19½"
Maximum Height.....	35½"
Width .....	27"
Length.....	76"



### HEAD REST ASSESSORIES

Crescent Removable Head Rest for Massage Therapy



### AM-SP220 MASSAGE THERAPY TABLE

The AM-SP220 Massage Therapy Table has 2-1/2" of multi layer foam padding. A firm density bottom layer with a soft top layer for patient comfort. The table is equipped to accept an adjustable crescent shaped headrest in addition to the standard headrest. The table is shown with Patina Vinyl, which has an ultra soft leather feel, superior wear resistance and a soft hi-loft backing.

Two piece top section .....	27" x 76".	Minimum Height.....	19½"
2½" Multi Layer Foam.		Maximum Height.....	35½"
Bottom Layer 1½" Firm Density Foam		Width .....	27"
Top Layer 1" Soft Density Foam.		Length.....	76"
One Piece Head Section 15" x 27" Raises 45 degrees. Lowers 90 degrees.			
One Piece Body Section 27" x 61".			

### AM-SP 225 ELEVATION / ADJUSTING TABLE

Four piece top section .....	72" x 24"
Minimum Height.....	19-1/2"
Maximum Height.....	35-1/2"
Width .....	24"
Length with footrest in .....	72"
Length with footrest extended ...	77"
Headrest .....	14" x 13"
Arm rest.....	19" x 5"
Face paper holder and cutter	

