



AM-SP200 TREATMENT TABLE

Two piece top section27" x 76"
 Minimum Height.....19½"
 Maximum Height.....35½"
 Width27"
 Length.....76"
 One piece head section.....15" x 27"
 Raises 45 degrees. Lowers 90 degrees.
 Contoured face/nose opening.
 One piece body section 61" x 27"



AM-SP250 TREATMENT TABLE

Four piece top section27" x 76"
 Minimum Height.....19½"
 Maximum Height.....35½"
 Width27"
 Length.....76"
 Three piece head section.....
 Center section.....15" x 13"
 Arm rest.....12½" x 5½"
 Raises 45 degrees. Lowers 90 degrees.
 Contoured face/nose opening.
 One piece body section61" x 27"



AM-SP350 TREATMENT TABLE

Three piece top section27" x 76"
 Minimum Height.....19½"
 Maximum Height.....35½"
 Width27"
 Length.....76"
 One piece head section.....15" x 27"
 Raises 45 degrees. Lowers 90 degrees
 Contoured face/nose opening.
 One piece center section22" x 27"
 Does not raise.
 One piece foot section.....39" x 27"
 Raises 70 degrees

