

**AM-SP400 TRACTION TABLE**



- Four piece top section .....27" x 76"
- Minimum Height.....20.5"
- Maximum Height.....36.5"
- Width .....27"
- Length.....76"
- Overall length .....95"
- One piece head section.....15" x 27"  
 Raises 45 degrees. Lowers 90 degrees.
- Contoured face/nose opening
- Thoracic section.....17 1/2" x 27"
- Small center section.....10" x 27"
- Foot section.....33" x 27"  
 Raises 80 degrees

Friction Free Lumbar Section separates 5" and rolls on steel ball roller bearings. The rolling Section can be locked in any position with a positive gear rack system.

**AM-SP450 TRACTION TABLE**



- Six piece top section .....27" x 76"
- Minimum Height.....20.5"
- Maximum Height.....36.5"
- Width .....27"
- Length.....76"
- Overall length .....95"
- Three piece head section
- Center section.....15" x 13"
- Arm rest.....12 1/2" x 5 1/2"  
 Raises 45 degrees. Lowers 90 degrees.
- Contoured face/nose opening.....17 1/2"
- Thoracic section.....17" x 27"
- Small center section.....10" x 27"
- Foot section.....33" x 27"  
 Raises 80 degrees

Friction Free Lumbar Section separates 5" and rolls on steel ball roller bearings. The rolling Section can be locked in any position with a positive gear rack system.

**AM-SP500 TREATMENT TABLE**



- Five piece top section .....27" x 76"
- Minimum Height.....19 1/2"
- Maximum Height.....35 1/2"
- Width .....27"
- Length.....76"
- Three piece head section .....
- Center section.....15" x 13"  
 Raises 45 degrees. Lowers 90 degrees.
- Contoured face/nose opening.
- One piece center section .....22" x 27"  
 Raises 30 degrees and has a spring loaded lock that must be released before the section will raise.
- One piece foot section.....39" x 27"  
 Raises 70 degrees