

AARON 950™-G OB/GYN TOTAL SYSTEM SOLUTION

We've combined our Aaron 950 package [see page 6] with our superb smoke evacuation system [see page 12] and our complete mobile stand [see page 22] to produce a compact, mobile system for the office or hospital. As an added benefit, we included a loop electrode kit (5 of our most popular loops) and a smoke evacuator hose reducer fitting [see page 13] for attaching to a vaginal speculum to create the most cost effective GYN system on the market today.

STYLE

A950G

STANDARD EQUIPMENT

- **A950 electrosurgical generator**
 - + Five extra disposable electrodes
 - + 50 sharp and 50 blunt disposable dermal tips
 - + A901 reusable single 3-button pencil
 - + A910-P ten handpiece drapes
 - + A1202C reusable grounding cord
 - + A1202 five solid grounding pads
- **SEGYN Smoke Shark system**
 - + SE01 Smoke Shark
 - + SF18 18-hour filter
 - + 786T 7/8" tube 6' long
 - + FSSE footswitch
 - + SERF reducer fitting
[attaches to vaginal speculum]
- **A812-C electrosurgical mobile stand**
 - + A812 mobile stand
 - + A812-E height extender
 - + A808-T instrument tray
 - + A812-BT bottom tray
- **ESLK electrosurgical loop kit**
 - + ES09 10 mm x 10 mm loop
 - + ES11 15 mm x 10 mm loop
 - + ES12 20 mm x 8 mm loop
 - + ES13 20 mm x 15 mm loop
 - + ES16 10 mm x 8 mm square

