

EPIC[®] MAGNIFIER

The perfect combination of world-class optics, precision movement and a large viewing area.

Burton's new Epic[®] Magnifier commands attention based on aesthetics alone – but its real value lies in its world class optics, performance and precision movement. With its patented arm system, the Epic[®] is more rugged than any other magnifier. The arm stops exactly where positioned and does not require constant tightening or lock mechanism adjustments.



PERFORMANCE RELIABILITY VALUE



Durable, precise, and powerful, Burton's Epic[®] Magnifier sets new performance standards in health care vision.

The Right Light

MAGNIFIERS

EPIC[®] MAGNIFIER

Precise, powerful, sturdy...and easy on the eyes.

EPIC[®] MAGNIFIER FEATURES:

- > 3.5 diopter polished glass lens
- > 6.5" (16.5 cm) diameter circular viewing area
- > Powerful 22-watt bulb – 10,000-hr. average bulb life
- > Parallel motion, patented enclosed arm system
- > Quick start electronic ballast (flicker-free)
- > Integral handles and cover for lens protection
- > Available in floorstand and table/wall models
- > 5-year warranty



MODELS

115V		230V	
#EM60FL	#EM62FL	Epic [®] Magnifier Floorstand on Mobile Base (casters included)	
#EM60W	#EM62W	Epic [®] Magnifier Wall/Table	

ACCESSORIES

#0009800PK	Set of 2 Replacement Blubs, 115V
#0009400PK	Set of 2 Replacement Blubs, 230V
#8970-001	4 Diopter Lens Attachment
#8970-002	6 Diopter Lens Attachment



Luxo Medical Products