

TX[®] Traction Tables

ADP-400 Traction Table

REF: ADP400

- Four sections
- Friction free gliding lumbar section
- Turret mounted traction pedestal allows for angular pulls
- Head and foot sections incline up to 90°
- Head section lowers to 25°
- 400 lb (181 kg) lifting capacity
- Adjustable height range between 21.5" to 41.5" (55 cm to 105 cm)
- Standard hand control
- Optional foot control
- Optional retractable caster
- Available in all upholstery colors
- Traction unit sold separately



Technical Specifications

Mains Power: 120 V, 50 Hz
Current: 1.6 A
Power Input: 200 W
Weight: 140 lbs (64 kg)
Shipping Weight: 200 lbs (91 kg)
Dimensions: 77" L x 27" W x 21.5" to 41.5" H
(196 cm x 69 cm x 55 cm to 105 cm)
Electrical Cass: Class 1, Type B



TTFT-200 Traction Table

REF: 6840

- Two sections
- Traction unit pedestal
- 350 lb (159 kg) capacity
- Fixed height of 32" (81 cm)
- Available in all upholstery colors
- Traction unit sold separately



TXS-1 Flexion Stool

REF: 7050

Adjustable height feature from 13"-20" (33-51 cm) allows therapist to place the patient's lumbar spine into neutral or flexed position for increased comfort. Upholstery can be color-coordinated to match traction table.

