

Ergostyle™



ErgoBasic™

REF: 9040 (soft foam), 9050 (firm foam)

- Tilt headpiece
- Accelerator™ Pelvic Drop
- Paper roll holder and cutter
- Height adjustable: 18"-24" (46-61 cm) – Ships at 22" (56 cm)
- Unitized steel frame
- Seamless cushions

Options:

- 89288** – E-Z Tilt headpiece
- 89295** – Gas spring-powered elevating pelvic section
- 89493** – Forward and toggle drop with tilt and elevation

Technical Specifications

Shipping Dimensions: 80" L x 30" W x 24" H (203 cm x 76 cm x 61 cm)
Shipping Weight: 125 lbs (56 kg)



ErgoBench™

- REF: 9000 – Tilt Headpiece (soft foam)
- REF: 9010 – Tilt Headpiece (firm foam)
- REF: 9020 – Fixed Top (soft foam)
- REF: 9030 – Fixed Top (firm foam)

Features

- Unitized Steel Frame
- Seamless Cushions
- Multi-fit Arm Rests
- Paper Roll Holder and Cutter
- Bottleneck Shape
- Height Adjustable: 18"-24" (46-61 cm) – Ships at 22" (56 cm)
- Available with Tilt headpiece



Technical Specifications

Shipping Dimensions: 80" L x 30" W x 24" H (203 cm x 76 cm x 61 cm)
Shipping Weight: 125 lbs (56 kg)

