

02 HEAT THERAPY

Mobile Heating Units

2402 – M-2 Mobile

33" H x 27" L x 16" W
(84 cm H x 69 cm L x 41 cm W)
Includes 12 Standard Size HotPacs

2502 – M-4 Mobile

33" H x 35" L x 20" W
(84 cm H x 89 cm L x 51 cm W)
Includes 24 Standard Size HotPacs

2302 – SS-2 Mobile

33" H x 21" L x 16" W
(84 cm H x 53 cm L x 41 cm W)
Includes 8 Standard Size HotPacs

Stationary Heating Units

2202 – SS

20" H x 18" L x 16" W
(51 cm H x 46 cm L x 41 cm W)
Includes 8 Standard Size HotPacs

2102 – E-1

16" H x 13" L x 8" W
(41 cm H x 33 cm L x 20 cm W)
Includes 4 Standard Size HotPacs

2105 – E-1 MHT

16" H x 13" L x 8" W
(41 cm H x 33 cm L x 20 cm W)
Includes 4 Hand Size HotPacs

2115K – E-1 MHT System (not shown)

16" H x 13" L x 8" W
(41 cm H x 33 cm L x 20 cm W)
Includes 4 Hand Size HotPacs,
2 Hand Contour Terry Covers and
2 Dual Pocket Hand Terry Covers

2802 – E-2

20" H x 15" L x 13" W
(51 cm H x 38 cm L x 33 cm W)
Includes 2 Oversize, 1 Neck
Contour and 3 Standard Size
HotPacs

Heating Unit Accessories

4010 – Side Table Rack

Fits Mobile Heating and Chilling
Units – 30" H x 20" L x 14" W
(76 cm H x 66 cm L x 35 cm W)
(No racks available for units
manufactured before 11/93)

4230 – Extra Shelf for Side Table Rack (not shown)

4228 – Dial Thermometer (not shown)

4226 – Automatic Time Switch (not shown)

4016 – Wall Mounted Towel Rack Size: 2" x 33" (5 cm x 84 cm) Six hooks extending 13" (38 cm) from wall (not shown)

HYDROCOLLATOR® HEATING UNITS

Hydrocollator® Heating Units

Hydrocollator Heating Units remain the benchmark against which all others are judged. Durable and easy to maintain, these high quality stainless steel units have an almost indefinite life and ensure a constant supply of temperature consistent HotPacs.



Side Table Rack

Stainless steel work surface and two wire shelves add efficiency to units. Use with models M-4, M-2 and SS-2.



Hydrocollator® Heating Units include:

- ◆ High-quality stainless steel.
- ◆ Thermostatically controlled temperature.
- ◆ 3" (8 cm) swivel-type rubber casters for silent, friction-free movement of mobile units.
- ◆ Easy maintenance. Simple to fill and drain. No plumbing required.
- ◆ UL listed and cUL listed (120 Volt).

Model E-1 MHT

- ◆ Designed for hand therapists.
- ◆ Non-insulated.
- ◆ Fits on countertop or mobile stand.



E-1 MHT

PARA-CARE® PARAFFIN BATH

Para-Care® Paraffin Bath

Para-Care provides penetrating heat that relieves pain and stiffness, and is optimal for patients suffering from arthritis, bursitis, muscle spasms and chronic joint inflammation. Para-Care's low-profile design allows easy access for the hand, elbow or foot, making it ideal for both home and clinical use. Two year warranty. ETL and cETL listed

See next page for accessories and ordering information.

