

LAB X SERIES

FEATURES

- 1 1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Armrest locks in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding throughout
- Neutral Grey powder-coated frame
- 2 knobs lock chrome armrest into position
- Adjustable foot levelers
- Easy-clean durable antimicrobial vinyl
- Complete series of chairs to choose from
- 400 lbs. capacity



Padded, adjustable flip-arm



2 knobs lock chrome armrest into position

	Seat	Height
6020-F	17" x 16"	20"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm and Drawer

- Adjustable, padded and upholstered flip-arm
- Easy-clean, one-piece drawer





	Seat	Height
6010-F	17" x 16"	20"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	27"—33"	

Blood Drawing Chair with Flip Arm

- Adjustable, padded and upholstered flip-arm



	Seat	Height
6022-F	17" x 16"	26"
	Width	Depth
Overall	43"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip Arm and Drawer

- Adjustable, padded and upholstered flip-arm
- Easy-clean, one-piece drawer
- Extra height reduces bending for staff
- Easy for patient to sit and exit chair
- Patient footrest

	Seat	Height
6011-F	17" x 16"	26"
	Width	Depth
Overall	34"	26"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall Blood Drawing Chair with Flip Arm

- Adjustable, padded and upholstered flip-arm
- Extra height reduces bending for staff
- Easy for patient to sit and exit chair
- Patient footrest



Arm in up position



	Seat	Height
6000X-F	35" x 18"	20"
	Width	Depth
Overall	48"	28"
	Height Adjustment	
Arms	29"—34"	

One arm in up position

Extra-Wide Blood Drawing Chair with Flip-Arms

- Padded and upholstered seat, back and armrests
- 2 adjustable flip-arms
- Accommodates large patient or child with parent up to 400 lbs. combined weight



	Seat	Height
6099-F	35" x 18"	26"
	Width	Depth
Overall	48"	28"
	Height Adjustment	
Arms	34"—40"	

One arm in up position

Extra-Tall & Wide Blood Drawing Chair with Flip-Arms

- Padded and upholstered seat, back and armrests
- 2 adjustable flip-arms
- Extra height reduces bending for staff
- Accommodates large patient or child with parent up to 400 lbs. combined weight



	Seat	Height
6000X-FD	35" x 18"	20"
	Width	Depth
Overall	58"	28"
	Height Adjustment	
Arms	29"—34"	

Extra-Wide Blood Drawing Chair with Drawer and Flip-Arm

- Padded and upholstered seat, back and armrests
- 1 adjustable flip-arm
- Easy-clean, one-piece drawer
- Accommodates large patient or child with parent up to 400 lbs. combined weight



	Seat	Height
6099-FD	35" x 18"	26"
	Width	Depth
Overall	58"	28"
	Height Adjustment	
Arms	34"—40"	

Extra-Tall & Wide Blood Drawing Chair with Drawer and Flip-Arm

- Padded and upholstered seat, back and armrests
- 1 adjustable flip-arm
- Easy-clean, one-piece drawer
- Extra height reduces bending for staff
- Accommodates large patient or child with parent up to 400 lbs. combined weight