



# ENTHERMICS Medical Systems

## EC770l

### Fluid Warming Cabinet



#### Specifications:

The single-chambered warming cabinet is constructed with 20 gauge stainless steel exterior casing and door. Window door and interior lighting allow observation of inventory with the door closed and convenient illumination when the door is open. The cabinet is warmed using the patented Halo Heat® low-heat-density electrothermal cable array. The cabinet is equipped with three (3) white, epoxy-coated wire baskets, each with a 14 liter capacity to accommodate most commonly used surgical fluids packaged in bags or bottles. The cabinet is furnished with a full-perimeter rubber bumper assembly and one (1) set of 5" (127mm) heavy duty casters — two with locking brake.

#### Special Information:

- Right-hand door swing is standard.
- Left-hand door is available as a special order.
- Stacking combinations listed under Options and Accessories.
- When ordering a combination of units to stack:
  - 1.) Indicate top and bottom placement.
  - 2.) Order one set of stacking hardware.
  - 3.) Delete full perimeter bumper (44092) and casters (CS-24874 and CS-24875) from upper unit.

#### Control:

The chamber is controlled by one (1) electronic control which consists of a 4 digit L.E.D. display, on/off key, heat indicator L.E.D, increase and decrease keys, built-in lock-out feature and a series of prompt sequence indicators. Individual temperature controls for irrigation fluids (IRR) and injection fluids (INJ) allow for precise temperature settings and can be easily set to operate in Celsius or Fahrenheit. The IRR temperature may be adjusted from 90° to 150°F (32° to 66°C), and the INJ temperature can be adjusted from 90° to 110°F (32° to 43°C). An alarm will sound if temperature exceeds 10°F (6°C) over the set-point temperature. In the event of a power failure, the cabinet will remember its programming and begin to operate as before. A warming shut-off system, separate from the electronic control, is an additional safety feature.

#### Capacity and Weights

Interior Capacity:	42 liters
Net Weight:	248 lb (112 kg)
Ship Weight:	325 lb (147 kg)

#### Electrical

120 V.A.C. 50/60 Hz, 1ph 1100 Watts, 9.0 Amps		NEMA 5-15P 15A - 125V Plug Hospital Grade	
230 V.A.C. 50/60 Hz, 1ph 1000 Watts, 4.5 Amps		BS 1363 Plug (UK only)	

#### Options and Accessories

Cylinder Lock for Door Handle	LK-22567
Leg Kit, 6" (152mm)	44093
Stacking Hardware:	CONTACT FACTORY
EC230 or EC230l over EC770l	
EC340 or EC340l over EC770l	
EC770 over EC770l	
EC770l over EC770l	
Window Door Assembly, Left Hand	5001255

