



Heavy Duty Narcotics Box, Dual Door/Double Lock, 24 x 16 x 8 #2703

Specifications



Construction:

- Heavy-duty 18 gauge steel

Features include:

- Dual Door Construction
- True double locking protection with two high security pick-resistant tubular locks
- Locks are keyed separately with two keys per lock included
- Unique lock sequences – all cabinets are keyed differently
- Reversible Cabinet can be positioned to open to either left or right
- Stainless steel pinned door hinges prevent pin removal
- 2 Adjustable height shelves
- Predrilled holes for mounting on a wall, cabinet or other solid surface (wood stud mounting hardware included)
- Narcotics box is delivered ready to use
- Five year warranty against manufacturer defect
- Durable hard-baked powder coat finish

Thirteen available Color Choices:

- White (-WH), Beige (-BG), Red (-RD), Light Blue (-LB), Navy (-NY), Mauve (-MV), Teal (-TL), Yellow (-YL), Light Gray (-GY), Burgundy (-BU), Hammer Tone Gray (-HTGY), Hammer Tone Blue (-HTDB), and Hammer Tone Green (-HTGN)

Overall dimensions:

- 24" H x 16" W x 8" D (61 cm H x 40.7 cm W x 20.4 cm D)
- Net weight: 32 pounds (14.6 kg)
- Approx. shipping weight: 37 pounds (16.9 kg)