



## Painted Steel Cast Cart Standard Package #6030



### Specifications

#### Features include:

- Painted Steel Body
- Aluminum Drawers
- Vinyl-padded full wrap-around bumper on a steel base
- Five-inch premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Heavy-duty ball-bearing drawer slides, warranted for the life of the cart
- Stainless steel bucket holder with vinyl mat and pump ole bracket
- Aluminum towel bar
- 24-inch stainless steel ruler, mounted on the top of the cart
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty

#### Drawer configuration:

- 1 ea 3” top drawer with black vinyl mat
- ea 3” drawers with four adjustable dividers per drawer
- 2 ea 6” drawers with four adjustable dividers per drawer

#### Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

#### Finished cart dimensions:

- 43” H x 44” W x 21” D (109.2 cm x 111.8 cm x 53.3 cm)
- Cart net weight: 182 lbs (82.6 kg)
- Approx. shipping weight 222 lbs (101 kg)



## Painted Steel Cast Cart Deluxe Package #6035



### Specifications

#### Features include:

- Stainless Steel Folding Shelf,
- Stainless Steel Side Storage Cabinet,
- Swivel Tourniquet / pump Pole
- Scissors cable
- Painted Steel Cabinet and Base
- Aluminum Drawers
- Vinyl-padded full wrap-around bumper on a steel base
- Five-inch premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Heavy-duty ball-bearing drawer slides
- Stainless steel bucket holder with vinyl mat and pump pole bracket
- Aluminum towel bar
- 24-inch stainless steel ruler, mounted on the top of the cart
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty

#### Drawer configuration:

- 1 ea 3” top drawer with black vinyl mat
- 5 ea 3” drawers with four adjustable dividers per drawer
- 2 ea 6” drawers with four adjustable dividers per drawer

#### Twelve available color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

#### Finished cart dimensions:

- 43” H x 44” W x 21” D (109.2 cm x 111.8 cm x 53.3 cm)
- With tourniquet / pump pole: 76” H x 44” W x 21” D (193 cm x 111.8 cm x 53.3 cm)
- Cart net weight: 209 lbs (95 kg)
- Approx. Shipping weight: 250 lbs (113.7 kg)