



## 4 Drawer Short Emergency Cart Standard Package Model #6300



### Specifications

#### Features Include:

- Short Cabinet with stainless steel breakaway locking system
- Breakaway lock seals, pkg-50
- Vinyl-padded full wrap-around steel bumper
- Five-inch premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Plastic removable top
- Stainless steel top rail/push handle
- Ball-Bearing Drawer Slides
- Double-walled drawer faces with no rough edges
- Cart is delivered fully assembled
- Certified "Latex Free"
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

#### Drawer Configuration:

- 2 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)
- 1 x 12-inch (30.5 cm)

#### Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

#### Twelve Available Color Choices:

Red (-RD), Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

#### Assembled Cart Dimensions:

38.5" H x 32" W x 22" D (97.2 cm H x 55.9 cm D x 81.3 cm W)

Cart weight: 119 pounds (54 kg)

Approx. shipping weight: 154 pounds (70 kg)