



6 Drawer Tall Medical Workstation Standard Package #6450

Specifications

Construction:

- Steel

Features include:

- Tall cabinet with key locking system
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) premium rubber single-wheel caster set, full swivel – one with brake, one directional
- Plastic removable top
- Aluminum push handles (2)
- Ball-bearing drawer slides
- Double-walled drawer faces with no rough edges
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish



Drawer configuration:

- 4 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)
- 1 x 12-inch (30.5 cm)

Drawer dimensions:

- 23” W x 17” D (58.4 cm W x 43.2 cm D)

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Assembled cart dimensions:

- 41.5” H x 32” W x 22” D (105.4cm H x 55.9cm D x 81.2cm W)
- Cart weight: 145 pounds (65.8 kg)
- Approx. shipping weight: 180 pounds (81.9 kg)