



## Secure Storage Anesthesia Cart with Key Lock Standard Package #6610



### Specifications

#### Construction:

- All steel

#### Features include:

- Tall cabinet with key lock system
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) sealed ball-bearing premium rubber single-wheel casters, full swivel – one with brake, one directional
- Aluminum push handle/side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Double-walled drawer faces with no rough edges
- Cart is delivered with wheels installed; some accessories require user installation.
- Certified “Latex Free”
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

#### Drawer configuration:

- 6 x 3-inch (7.6 cm)
- 2 x 6-inch (15.2 cm)

#### Drawer dimensions:

- 23” W x 17” D (58.4 cm W x 43.2 cm D)

#### Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

#### Assembled cart dimensions:

- 41.5” H x 32” W x 22” D (105 cm H x 81.3 cm W x 55.9 cm D)
- Cart weight: 120 pounds (54.4 kg)
- Shipping weight: 155 pounds (61.2 kg)



## Secure Storage Anesthesia Cart Specialty Package with Plastic Top #6611



### Specifications

#### Construction:

- All steel

#### Features include:

- Stainless steel top rail/push handle
- Utility hooks, set-4
- 6-Compartment drawer trays, quantity of 3
- 7-Compartment drawer trays, quantity of 3
- Tall cabinet with key locking system
- Vinyl-padded full wrap-around bumper
- Five-inch (12.7 cm) sealed ball-bearing premium rubber single-wheel casters, full swivel – one with brake, one directional
- Plastic removable top
- Ball-bearing drawer slides
- Double-walled drawer faces with no rough edges
- Cart is delivered with wheels installed; some accessories require user installation.
- Certified “Latex Free”
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

#### Drawer configuration:

- 6 x 3-inch (7.6 cm)
- 2 x 6-inch (15.2 cm)

#### Drawer dimensions:

- 23” W x 17” D (58.4 cm x 43.2 cm)

#### Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

#### Assembled cart dimensions:

- 44.5” H x 32” W x 22” D (113 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 125 pounds (56.7 kg)
- Shipping weight: 160 pounds (63.5 kg)