



Malignant Hyperthermia Cart Standard Package #6640



Specifications

Construction:

- Steel Cabinet and Drawers with Built-in Refrigerator

Features:

- 1.6 cu ft refrigerator/freezer with hospital grade plug
- 10 Glycol packs
- LED-illuminated digital thermometer with cart power safety feature.
- Separate breakaway locks on both the refrigerator and the two storage drawers
- Tall cabinet
- Breakaway lock seals, unnumbered, package of 50
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) sealed ball-bearing premium rubber single-wheel casters, full swivel – one with brake, one directional
- Aluminum push handle/side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty
- Compatible with many Harloff drawer dividers and accessories
- Hard-baked powder coat finish

Drawer configuration:

- 1 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)

Drawer dimensions:

- 23” W x 17” D (58.4 cm W x 43.2 cm D)

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Assembled cart dimensions:

- 41.5” H x 32” W x 22” D (105 cm H x 81.3 cm W x 55.9 cm D)
- Cart net weight: 146 pounds (66.4 kg)
- Approx. Shipping weight: 181 pounds (82.3 kg)



Malignant Hyperthermia Cart Specialty Package #6641

Specifications



Construction:

- Steel Cabinet and Drawers with Built-in Refrigerator

Features:

- 1.6 cu ft refrigerator/freezer with hospital grade plug
- 10 Glycol packs
- LED-illuminated digital thermometer with cart power safety feature.
- Separate breakaway locks on both the refrigerator and the two storage drawers
- IV Pole
- Stainless Steel Top Rail/Push Handle
- Side-mounted drop shelf
- Full Drawer Tray organizer with adjustable dividers for top drawer
- Tall cabinet
- Breakaway lock seals, unnumbered, package of 50
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) sealed ball-bearing premium rubber single-wheel casters, full swivel – two with brake
- Aluminum push handle/side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty
- Compatible with many Harloff drawer dividers and accessories
- Hard-baked powder coat finish

Drawer configuration:

- 1 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Drawer dimensions:

- 23” W x 17” D (58.4 cm W x 43.2 cm D)

Assembled cart dimensions:

- 44.5” H x 32” W x 22” D (113.0 cm H x 81.3 cm W x 55.9 cm D)
- IV Pole extends to 85” (216 cm)
- Cart weight: 154 pounds (70 kg)
- Approx. Shipping weight: 199 pounds (90.5 kg)