



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Classic Line, Six Drawer Workstation (Tall) with
 Overhead Shelving, Key Lock, Specialty Package
 #6650**



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Tall cabinet with key locking system – two keys
- Short overhead shelf system with two 14-inch deep (35.6 cm) overhead monitor shelves
- Extended width full wrap-around steel bumper for stability when shelves are loaded
- Extended width full wrap-around vinyl padded bumper increases stability when shelves are loaded
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Aluminum push handle / side rail
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered with wheels installed; some accessories require user installation.
- Full range of accessories and drawer organizers available
- Pre-drilled threaded holes for side mounted accessories
- Twelve year warranty
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (3) – 3" (7.6 cm)
- (2) – 6" (15.2 cm)
- (1) – 9" (22.9 cm)

Sixteen available color choices:

***Antimicrobial Finish Available**

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Internal Drawer Dimensions:

- 20" W x 17.25" D (50.8 cm W x 13.8 cm D)

Assembled Cart Dimensions

- 41.5" H x 32" W x 22" D (105.4 cm H x 81.3 W x 55.9 D)
- With Shelves: 66.5" (168.9 cm) H
- Cart weight: 180 pounds (81.7 kg)
- Shipping weight: 220 pounds (100 kg)

<http://www.ariamedical.com/harloff-6650.html>

Toll free: 800-330-3591
 Local: (210) 281-9602
 Fax: (210) 281-9604

email: sales@ariamedical.com
 website: <http://www.ariamedical.com>