



HARLOFF OptimAL™ LINE ALL-ALUMINUM PUNCH CARD MEDICATION CART #AL460PC



Left Column Drawers, 20.25" W x 17" D:

- 1 ea – 4" Multipurpose Storage Drawer
- 3 ea – 9.75" Punch card drawers with 3 adjustable row dividers

Right Column Drawers, 6.625" W x 17" D:

- 1 ea – 4" Top Multipurpose Storage Drawer
- 1 ea – 9.75" Locking Narcotics Drawer
- 2 ea – 9.75" Multipurpose Storage Drawer

*Using the 3 ea. - 9.75" drawers for punch cards allows an increase capacity of 120 Punch Cards

Assembled Cart Dimensions:

45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)

Cart Net weight: 137 lbs (62.2 kg)

Approx. Shipping weight: 197 lbs (89.6 kg)

Twelve available End Panel/Drawer Post color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Gray (-TXGY) and Hammertone Green (-TXGN)

Specifications

480 Punch Card Capacity (360 left column + 120 right column)

AL460PC KEY FEATURES AND BENEFITS:

- Spacious Medication Delivery System with 480 Punch Card Capacity
- The top of the cart is a smooth and sturdy surface for carrying regulatory documents and pamphlets
- Pull out writing shelf for additional writing space for logging patient's charts
- Easy to push All Aluminum Construction, with Integrated push handles for easy transporting
- Maximum Security, built with BEST® Brand Locks on both cabinet and narcotics drawer
- Ball-bearing, full extension, self-closing drawer slides
- The "slam action" lock on the Narcotics drawer allows the drawer to be closed without using the key. An internal drawer cover keeps the drawer contents secure, even if all other drawers are removed.
- Full wrap around vinyl bumper, protecting the cart
- 5" Sealed ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional, easy to maneuver over all floor types; tile, carpet, linoleum, wood and concrete
- Offered in a "Specialty Package," which includes drawer tray with dividers, for top miscellaneous drawer, waste container and locking sharps container, # AL460PCSP
- Hard-baked powder coat finish.



Optima^{AL}™ LINE ALL-ALUMINUM PUNCH CARD MEDICATION CART #AL600PC



Specifications

600 Punch Card Capacity (480 left column + 120 right column)

AL600PC KEY FEATURES AND BENEFITS:

- Spacious Medication Delivery System with 600 Punch Card Capacity
- The top of the cart is a smooth and sturdy surface for carrying regulatory documents and pamphlets
- Pull out writing shelf for additional writing space for logging patient's charts
- Easy to push All Aluminum Construction, with Integrated push handles for easy transporting
- Maximum Security, built with BEST® Brand Locks on both cabinet and narcotics drawer
- Ball-bearing, full extension, self-closing drawer slides
- The "slam action" lock on the Narcotics drawer allows the drawer to be closed without using the key. An internal drawer cover keeps the drawer contents secure, even if all other drawers are removed.
- Full wrap around vinyl bumper, protecting the cart
- 5" Sealed ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional, easy to maneuver over all floor types; tile, carpet, linoleum, wood and concrete
- Offered in a "Specialty Package," which includes drawer tray with dividers, for top miscellaneous drawer, waste container and locking sharps container, # AL600PCSP
- Hard-baked powder coat finish

Left Column Drawers, 27.5" W x 17" D:

- 1 ea – 4" Multipurpose Storage Drawer
- 3 ea – 9.75" Punch card drawers with 3 adjustable row dividers

Right Column Drawers, 6.625 W x 17" D:

- 1 ea – 4" Top Multipurpose Storage Drawer
- 1 ea – 9.75" Locking Narcotics Drawer
- 2 ea – 9.75" Multipurpose Storage Drawer
- *Using the 3 ea. - 9.75" drawers for punch cards allows an increase capacity of 120 Punch Cards

Twelve available End Panel/Drawer Post color choices:

Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Gray (-TXGY) and Hammertone Green (-TXGN)

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D (116.2 cm x 109.2 cm x 59.7 cm)
- Cart weight: 142 lbs (64.5 kg)
- Approx. shipping weight: 197 lbs (89.6 kg)