



Six Drawer MR-Safe Anesthesia Cart #MR6K



Specifications

Construction

- All-Aluminum Cabinet and Drawers

Standard Features

- Tall Cabinet (33.25 inches / 84.5 cm) with BEST® Lock keylocking system
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) MR-Safe single-wheel caster set, full swivel – two with brake
- Stainless Steel Pull-out Shelf
- Pre-threaded accessory mounting holes
- Low Ferrous Ball-Bearing Drawer Slides
- Replaceable Plastic Top with molded-in Push Handles
- Cart is delivered fully assembled
- Certified “Latex Free”
- Certified MR-Safe by independent testing authority (copies of report available)
- Five Year Warranty
- MR-Safe accessories available
- Hard-baked powder coat finish

Twelve available drawer color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Gray (-TXGY) and Hammertone Green (-TXGN)
- Shelf Painted Gloss White with MRI Sticker

Drawer configuration:

- 1 x 4-inch (10.2 cm)
- 2 x 3.25-inch (8.3 cm)
- 2 x 6.5-inch (16.5 cm)
- 1 x 9.75-inch (24.8 cm)

Assembled Cart Dimensions:

- 45.7 5” H x 30.5” W x 23.75” D (116.2 cm H x 77.5 cm W x 60.3 cm D)
- Cart Net Weight: 126 pounds (57.3 kg)
- Approx. Shipping Weight: 156 pounds (71 kg)



Six Drawer MR-Safe Anesthesia Cart Specialty Package #MR6K-MAN



Specifications

Construction

- All-Aluminum Cabinet and Drawers

Standard Features

- 3 Tiered Raised back rail system
- Side Rails (4)
- Full drawer tray organizers (2)
- Tall Cabinet (33.25 inches / 84.5 cm) with BEST® Lock keylocking system
- Vinyl-padded full wrap-around steel bumper
- Five-inch (12.7 cm) MR-Safe single-wheel caster set, full swivel – two with brake
- Stainless Steel Pull-out Shelf
- Pre-threaded accessory mounting holes
- Low Ferrous Ball-Bearing Drawer Slides
- Replaceable Plastic Top with molded-in Push Handles
- Cart is delivered fully assembled
- Certified “Latex Free”
- Certified MR-Safe by independent testing authority (copies of report available)
- Five Year Warranty
- MR-Safe accessories available
- Hard-baked powder coat finish

Twelve available drawer color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Gray (-TXGY) and Hammertone Green (-TXGN)
- Shelf Painted Gloss White with MRI Sticker

Drawer configuration:

- 1 x 4-inch (10.2 cm)
- 2 x 3.25-inch (8.3 cm)
- 2 x 6.5-inch (16.5 cm)
- 1 x 9.75-inch (24.8 cm)

Assembled Cart Dimensions:

- 45.75” H x 33” W x 25” D (116.2 cm H x 81.3 cm W x 63.5 cm D) 68.5” H (174 cm) with Back Rail System
- Cart Net Weight: 139 pounds (63.1 kg)
- Approx. Shipping Weight: 169 pounds (77 kg)