



## Standard Line 210 Capacity Unit Dose Medication Box Cart Standard Package #SL210BOX

### Specifications



#### Construction:

- Steel Body, Aluminum drawers

#### Features:

- Dual Column Medication Delivery Cart with 210 Medication Box capacity
- Dual Locking System complies with regulatory security standards
- BEST® Brand Lock on cabinet; Illinois Lock on narcotics drawer
- The “slam action” lock on the Narcotics drawer allows the drawer to be closed without using the key.
- An internal drawer cover keeps the narcotics drawer contents secure, even if all other drawers are removed.
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5” Sealed ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional, easy to maneuver over all floor types; tile, carpet, linoleum, wood and concrete
- Hard-baked powder coat finish
- Offered in a **Specialty Package** which includes pullout side shelf, drawer tray with dividers, for top miscellaneous drawer, waste container and locking sharps container - see **SL210BOXSP**

#### Left Column Drawers:

- 1 ea – 4” Top Multipurpose Storage Drawer with internal divider set, 30 box capacity per drawer
- 6 ea – 3.25” Med Box Drawers with internal divider set, 30 box capacity per drawer
- 1 ea – 9.75” Multipurpose Storage Drawer

#### Thirteen available color choices for side panels:

- Hammer Beige (HB), Gloss Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Navy (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)
- Standard Drawer color is Hammer Beige

#### Right Column Drawers:

- 1 ea – 4” Top Multipurpose Storage Drawer
- 1 ea – 9.75” Locking Narcotics Drawer
- 2 ea – 9.75” Multipurpose Storage Drawer

#### Assembled Cart Dimensions:

- 45.75” H x 36” W x 23.5” D (116.2 cm x 91.4 cm x 59.7 cm)
- Cart net weight: 198 lbs (90 kg)
- Approx shipping weight: 243 lbs (110.5 kg)