



Value Line Punchcard Medication Cart with Key Lock #VL280PC



Specifications

Construction:

- Steel

Features:

- Tall Cabinet with key locking system
- Pull out Side Shelf
- Capacity of up to 280 punch cards
- Bumper Frame with projected front casters and (4) Corner Bumpers
- Five-inch premium rubber single-wheel casters, full swivel – one with brake, one directional
- Painted steel top rail/push handle
- Removable plastic top
- Full extension Ball-Bearing Drawer Slides
- Cart is delivered fully assembled
- Certified “Latex Free”
- Five Year Warranty
- Compatible with most accessories
- Hard-baked powder coat finish

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Gray (-TXGY) and Hammertone Green (-TXGN)

Drawer configuration:

- 1 x 3-inch (7.6 cm) multipurpose storage drawer
- 2 x 9-inch (22.9 cm) punchcard drawers with two adjustable dividers
- 1 x 9-inch (22.9 cm) drawer with locking narcotics compartment and one adjustable divider

Drawer Dimensions:

- 23” W x 17” D (58.5 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 44.75” H x 29.5” W x 22.5” D (113.7 cm H x 75 cm W x 57.2 cm D)
- Cart weight: 140 pounds (63.7 kg)
- Approx. shipping weight: 175 pounds (79.6 kg)