


Safety bed rails are some of the most important security and accident-prevention accessories used in the patient room today. That's why all Invacare safety bed rails feature durable chrome-plated, welded steel construction, and are designed to withstand heavy-duty use. These rails can be installed permanently or semi-permanently on all spring-fabric medical-style beds, and they are easy to attach—no tools are necessary for assembly on most models.



Invacare® Economy Spring-Loaded Telescoping Bed Rail

Model no. 6650

- Durable 7/8" chrome-plated, welded steel construction
- Self-adjusting spring-loaded crossbraces snap easily into place without tools for quick installation
- **NEW** Third-location height adjustment permits middle positioning to eliminate the gap created between the mattress and rail when the head of the bed is raised (caregiver attendance still required)



Invacare® Spring-Loaded Telescoping Bed Rail

Model no. 6629


- Durable 7/8" chrome-plated, welded steel construction
- Self-adjusting spring-loaded crossbraces snap easily into place without tools for quick installation



Invacare® Deluxe Spring-Loaded Telescoping Bed Rail

Model no. 6628

- Heavy-duty, 1" chrome-plated, welded steel construction
- Self-adjusting spring-loaded crossbraces snap easily into place without tools for quick installation
- **NEW** Third-location height adjustment permits middle positioning to eliminate the gap created between the mattress and rail when the head of the bed is raised (caregiver attendance still required)



Invacare® Telescoping Home Bed Rail

Model no. 6631


- Tubular 7/8" chrome-plated, welded steel side rails have crossbraces with ten preselected lengths from 36" to 62"
- For use with home or dormitory-style beds; fits between mattress and box spring
- **NEW** Rails fit most single, double and queen-size beds (do not fit king-size beds)



Invacare® Half-Length Bed Rail

Model no. 6640

- Heavy-duty 1" chrome-plated, welded steel construction
- Single flat crossbrace is secured with bolt-on clamp assembly (assembly required)
- For optimal security and convenience, four half-length rails are recommended (one set for head area, one set for foot area)
- Use with home care or bariatric bed



Invacare® Half-Length Bed Rail

Model no. 6630

- Heavy-duty, 1" chrome-plated, welded steel construction
- Teflon®-coated bushings ensure smooth raising and lowering action
- Dual flat-channel crossbraces are secured with wing-nut clamp assemblies (some assembly required)
- For optimal security and convenience, four half-length rails are recommended (one set for head area, one set for foot area)

Specifications

Model no.	Interchanges with Model nos.	Rail length		Rail height above deck		Deck-width adjustability		Width between bed rails	Shipping weight	Warranty
		min.	max.	raised	lowered	min.	max.			
6623	6627 & 6629	47"	74"	14.5"	1.75"	32"	36"	39"	29 lb.	3 years limited
6625	6658 & 6628	47"	74"	14.5"	.50"	32"	36"	39"	35 lb.	3 years limited
6628	6625 & 6658	47"	74"	15.25"	1"	31"	37"	34.5" - 40.5"	35 lb.	3 years limited
6629	6623 & 6627	47"	74"	15.25"	1"	31"	37"	34.5" - 40.5"	29 lb.	3 years limited
6630		30"	30"	15"	2"	32"	36"	38"	29 lb.	3 years limited
6631		47"	74"	16.25"	3.75"	36"	62"	36" - 62"	29 lb.	3 years limited
6640		30"	30"	15"	1.50"	34"	36"	37" - 39"	25 lb.	2 years limited
6650		47"	74"	15.25"	1"	31"	37"	35" - 41"	29 lb.	2 years limited

NOT SHOWN: Model no. 6623 Invacare® Clamp-On Telescoping Bed Rail, Model no. 6625 Invacare® Deluxe Clamp-On Telescoping Bed Rail.

Invacare overbed tables are practical, easy-to-use units suitable for both home and institution. The table height easily adjusts for comfort while reading, writing or eating. These tables feature an attractive walnut wood-grained laminate top with a T-molded edge that is designed to resist damage. The heavy gauge chrome-plated, tubular welded steel construction ensures durability, and the 1.5" casters ensure smooth mobility of the table in any direction.



Invacare® Auto-Touch Overbed Table

Model no. 6417

- Light upward touch on patented spring-loaded locking handle allows table to elevate freely
- Light downward touch on table top locks table in position
- Unique "U" base design permits use with chairs



Invacare® Tilt-Top Overbed Table

Model no. 6418

- Adaptable tilt-top allows three angle positions in either direction using tilt-release lever
- Light upward touch allows table to elevate freely
- Light downward touch on table top locks it into position



Invacare® CareGuard™ Overbed Table

Model no. 6400

- Light upward touch allows table to elevate freely
- Light downward touch on table top locks table in position
- Economically priced

Specifications

Model no.	6417	6418	6400
Dimensions			
Table top	.75" h x 15" w x 30" l	.75" h x 15" w x 30" l	.75" h x 15" w x 30" l
Base	15.5" w x 26" l	15.5" w x 26" l	16.5" x 26.25"
Height adjustment	29" - 45"	25.75" - 39"	27" - 42"
Weight	27 lb.	28 lb.	22 lb.
Column/base/support bracket	14/16/13 gauge	14/16/13 gauge	14/14/13 gauge
Warranty	3 years limited	3 years limited	2 years limited