

Invacare trapeze units are an important patient room accessory designed to help patients change positions while in bed and to aid in transfer from bed to chair with minimum attendant assistance. Trapeze units have a wide range of height adjustments and handbar positions to maximize patient accessibility. The unit can be installed quickly and easily on an Invacare floor stand or on any Invacare home care bed end.



7714 Shown with #7740 Trapeze

### Invacare® Trapeze Floor Stand

Model no. 7714

- Provides free-standing support for Model no. 7740 fixed trapeze unit
- Folds flat for convenient storage and handling
- Assisted weight capacity of 250 lb.



### Invacare® Fixed Offset Trapeze Bar


Model no. 7740

- Heavy-gauge, octagon steel construction; no tools required for easy assembly
- Wall bumpers and padded mounting brackets protect bed ends and wall surfaces
- For use on Model no. 7714 floor stand or on metal frame headboards only
- Assisted weight capacity of 250 lb.
- Includes bar, brackets and grab bar

### Specifications

Model no.	7740	7714
Dimensions		
Overall height	61"	38"
Octagon steel tubing, diameter	1.25"	1.25"
Reach	34"	n/a
Base width	n/a	32"
Base length	n/a	40"
Weight	17 lb.	24 lb.
Shipping	via truck	via UPS
Warranty	3 years	3 years

## Bed Accessories



### Bed Buggy

Model no. 5140

- Used to transport split spring or assembled bedspring quickly and easily
- Allows for one-person set-up in much less time
- Can be used on all Invacare beds



### Extender

Model no. 5143

- Designed to expand 80" bed frame to 84"
- Complete with extended head pull tube and drive shaft, link fabric and helicals
- Not available on bariatric bed Model no. BAR600



### Stor-A-Bed

Model no. 5142

- The ultimate bed storage unit. Uses both bed ends to transport or store bed. Helps maintain bed inventory, requiring less warehouse space
- Can be used on all Invacare beds