

The Invacare® trapeze unit is an important patient room accessory designed to help patients change positions while in bed and to aid in transferring from bed to chair with minimum attendant assistance. Trapeze units have a wide range of height adjustments and handbar positions to maximize patient accessibility.



7714P Shown with #7740P Trapeze

Invacare® Trapeze Floor Stand

Model no. 7714P

- Provides free-standing support for model no. 7740P fixed trapeze unit
- Folds flat for convenient storage and handling
- Assisted weight capacity of 250 lb.



Invacare® Fixed Offset Trapeze Bar

Model no. 7740P

- Heavy-gauge, octagon steel construction; no tools required for easy assembly
- Wall bumpers and padded mounting brackets protect bed ends and wall surfaces
- For use on model no. 7714P floor stand or on any Invacare home care bed end
- Assisted weight capacity of 250 lb.
- Includes bar, brackets and grab bar

Specifications

Model no.	7740P	7714P
Dimensions		
Overall height	61"	38"
Octagon steel tubing, diameter	1.25"	1.25"
Reach	34"	n/a
Base width	n/a	32"
Base length	n/a	40"
Product weight	17 lb.	24 lb.
Shipping	via truck	via UPS
Limited warranty	3 years	3 years

Invacare® Bariatric Trapeze

Model no. BARTRAP

- Three locking overhead boom angle positions for comfortable, secure use
- Tool-free design simplifies assembly
- Four piece sub-assembly significantly reduces setup time
- Lightweight design for easy delivery
- Free-standing structure permits use with Bariatric, Home Care or non-medical bed
- 1000 lb. client weight capacity

Specifications

Model no.	BARTRAP
Dimensions	
Overall height	73"
Reach	35.5"
Base width	40"
Base length	50.5"
Product weight	72 lb.
Shipping weight	80 lb.
Limited warranty	3 years

Invacare, the Blue Medallion Design and Yes, you can are registered trademarks and Bariatric Solutions Design, IVC and CareGuard are trademarks of Invacare Corporation. Specifications are subject to change without notice.

© 2005

Invacare Corporation
Form No. 00-274
Rev. 10/05