

HTR Tilt and Recline Chairs

Invacare's HTR tilt and recline chairs provide your patients with exceptional comfort, positioning and support. Each model in the HTR line of chairs is built to last, and is easy to adjust and maintain. They can also be configured with a number of accessories, making them an excellent value for any long-term care facility that requires optimal positioning for patient care.

HTR Standard Features*

- Tilt adjustment (0° - 35°) for greater positioning
- Recline adjustment (95° to 125°) for improved positioning
- Fold-down back facilitates transfers and personal grooming
- Sliding seat lock for easy seat-depth adjustment (16" - 20")
- Flip-back adjustable height arms (7.5" - 12.5")
- Available in frame widths of 16" or 18"
- Available in three colors: mauve, grotto and blueberry

HTR 3000, HTR 5000 Options

- Three styles of front riggings (Model Nos. HRA92, HRA93, HRA95)

HTR 5000, HTR 5500 Options

- 20" seat-to-floor height option with 20" rear wheel and 6" front caster
- Width reduction of 4" utilizing 12" rear wheels with flat free inserts

* See back page for HTR specifications

Available Upholstery Colors



#427 Blueberry



#425 Mauve



#426 Grotto



HTR 3000 Features

Model no. HTR3000

- Three swivel-locking and one directional-locking 5" casters
- Sliding seat pan for use with pressure-relief cushions



HTR 3500 Features

Model no. HTR3500

- Three swivel-locking and one directional-locking 5" casters
- Contoured seat cushion
- Head wings
- Contoured headrest
- Lap tray
- One-piece elevating footboard



HTR 5000 Features

Model no. HTR5000

- 24" rear wheels and 8" casters
- Seat pan for use with pressure-relief cushions



HTR 5500 Features

Model no. HTR5500

- 24" rear wheels and 8" casters
- Contoured seat cushion
- Head wings
- Contoured headrest
- Lap tray
- One-piece elevating footboard

HTR Accessories

Substitute XX in part numbers below with number corresponding to your color choice: 25-Mauve, 26-Grotto and 27-Blueberry. Part numbers for the accessories include the hardware kits to attach to the chairs



Headrest Extension
(Option No. HRA1)
(Part No. 460069XX)
Dimensions 6" x 15"
Adjustment Range 3"
Hardware included



4" Abduction Pommel
(Option No. HRA2)
(Part No. 460070XX)
Hardware included
Not available on
HTR3000/5000




2" Side Bolsters (Pair)
Option No. HRA10
Part No. 450148XX
Dimensions 6" x 15"
Reduces seat width



Lateral Trunk Support (Pair)
(Option No. HRA20)
(Part No.:
1072292 – Mauve
1072299 – Grotto
1072300 – Blueberry)
Dimensions 4.5" x 6"
Adjustment Range 4"
Hardware included



Seat Positioning Strap
– Auto Style
(Option No. HRA30)
(Part No. 450149)
Hardware included



Adjustable Height IV Pole
(Option No. HRA35)
(Part No. 1072293)
Adjustment Range 35" to 61"
Hardware included



Standard Lap Tray
(Option No. HRA42)
(Part No. 460087 - 16")
(Part No. 460088 - 18")
Dimensions:
16" wide - 23" x 18"
18" wide - 23" x 20"
Standard on HTR3500
and HTR5500



Vinyl Tray Pad 1"
(Option No. HRA4CA)
(Part No. 460071XX - 16")
(Part No. 460072XX - 18")
Secures with Hook and
Loop Straps



Contoured Seat
Option No. HRA43
Part No. 450088-XX - 16"
Part No. 450089-XX - 18"
Standard on HTR3500 and
HTR5500




Headwings (Pair)
(Option No. HRA44)
(Part No.:
1072294 – Mauve
1072301 – Grotto
1072302 – Blueberry)
Dimensions 7" x 11.5"
Adjustment Range 4"
Hardware included
Standard on HTR3500
and HTR5500



Winged Headrest
(Option No. HRA46)
(Part No. 450127-XX)
Secures with hook and
loop straps



Contoured Headrest
(Option No. HRA47)
(Part No. 450094-XX - 16")
(Part No. 450095-XX - 18")
Secures with hook and
loop straps



Front Anti-Tippers (pair)
(Option No. HRA60)
(Part No. 460073)
Add 9" to chair depth



Footboard Pad
(Option No.:
HRA68-16 - 16")
HRA68-18 - 18")
(Part No. 450121 - 16")
(Part No. 450122 - 18")
Secures with hook and
loop straps



12" Wheel with Flat
Free Inserts
(Option No. HRA12)
Optional on HTR5000
and HTR5500
Width reduction of 4"



Swing Away Footrests
(Option No. HRA92)
5.5" Horizontal adjustment
5" Footboard adjustment
15" Total additional depth



Articulating Elevating
Leg Rest
(Option No. HRA95)
5.5" Horizontal adjustment
13" Total footboard
adjustment
29" Total additional depth



One-Piece Elevating Footboard
(Option No. HRA93)
5.5" Horizontal adjustment
7.5" Footboard adjustment
23" Total additional depth
Standard on HTR3500
and HTR5500

HTR Specifications

Model Nos. HTR3000, HTR3500, HTR5000, HTR5500

	HTR3000	HTR3500	HTR5000	HTR5500
Overall height	50"	50"	50"	50"
Overall width	22"/24"	22"/24"	26"/28"	26"/28"
Overall depth upright	39"	57"	39"	57"
Overall depth recline	52"	69"	52"	69"
Seat width	16" or 18"	16" or 18"	16" or 18"	16" or 18"
Seat depth	16"-20"	16"-20"	16"-20"	16"-20"
Seat-to-back height	34"	31"	34"	31"
Seat-to-floor height	22"	22"	20" or 22"	20" or 22"
Tilt range	0 - 35°	0 - 35°	0 - 35°	0 - 35°
Recline range	95 - 125°	95 - 125°	95 - 125°	95 - 125°
Overall tray (WxL)	16" or 18" x 23.25"	16" or 18" x 23.25"	16" or 18" x 23.25"	16" or 18" x 23.25"
Weight capacity	250 lb.	250 lb.	250 lb.	250 lb.
Product weight (16/18 wide)	66 lbs./69 lbs.	98 lbs./101 lbs.	88 lbs./91 lbs.	105 lbs./108 lbs.
Shipping weight (16/18 wide)	114 lbs./117 lbs.	146 lbs./149 lbs.	136 lbs./139 lbs.	153 lbs./156 lbs.
Carton dimensions (LxWxH)	46" x 33" x 41.5"	46" x 33" x 41.5"	46" x 33" x 41.5"	46" x 33" x 41.5"
Warranty on chair	3years	3years	3years	3years
Warranty on upholstery	1 year	1 year	1 year	1 year



HTR 5500 in full tilt position



HTR 5500 in full recline position



HTR 5500 in full tilt and recline position