

Invacare®
Bariatric Trapeze
and Related
Bariatric Bed
Products

Invacare® Bariatric Trapeze 
and Related Bariatric Bed Products



Yes, you can.™



Yes, you can.™

Invacare's new Bariatric Trapeze is designed to assist individuals weighing up to 600 pounds in changing positions while in bed and aid in transferring in and out of bed. In addition to providing maximum support and comfort, this innovative lightweight, tool-free design makes delivery and assembly simple. With its ease of operation and delivery, the Invacare® Bariatric Trapeze offers comfort to both clients and providers alike.



Model no. BARTRAP

Features

- 600 lb. client weight capacity
- Three locking overhead boom angle positions for comfortable, secure use
- Tool-free design simplifies assembly
- Four piece sub-assembly significantly reduces set-up time
- Lightweight design for easy delivery
- Free standing structure permits use with Bariatric, Home Care or Non-medical Bed

Specifications


Model no. BARTRAP

Dimensions	
Overall height	73"
Reach	35.5"
Base width	40"
Base length	50.5"
Product weight	72 lb.
Shipping weight	80 lb. / UPS
Client weight capacity	600 lb.
Limited warranty	3 years

Invacare is a registered trademark and Yes, you can. is a trademark of Invacare Corporation. Specifications are subject to change without notification.

Invacare Bariatric Bed

Invacare's Bariatric Bed is an extra wide, heavy-duty design capable of supporting up to 600 pound client weight. This full electric bed provides convenient motorized operation of the back and knees, as well as the bed height. Even with the heavy-duty construction, the BAR600 split-spring design makes delivery and storage easy. Invacare's Bariatric Mattress and Rails are designed to complement the bed for comfortable client use. Designed to meet KO549.



A Invacare Bariatric Bed (model no. BAR600), **B** Bariatric Mattress (model no. BARMATT) and **C** Bed Rails (model no. 6640)