


The Invacare Group I Therapeutic Support Surfaces are designed to reduce patient interface pressure for high risk individuals to aid in the treatment of pressure ulcers. Four types of modalities allow you to choose the right product for the individual user. Product types include alternating pressure, gel or static air overlays, and foam mattress replacements. All of these products are easy to deliver, set up and simple to use.



Invacare® CareGuard™ SA Static Air Overlay

Model no. CGSA1100

- 18 gauge, anti-microbial, medical grade vinyl
- Three-inch air cell height to reduce patient interface pressure
- Two inflation valves for use with hand pump (included) or electric inflator
- CPR release valve for rapid deflation required in the ICU
- Easy-to-adjust straps secure the overlay to any size mattress
- Designed to meet EO197
- Does not contain latex



Invacare® CareGuard™ 101 Home Care Foam Mattress

Model no. CG10180

- Five zones allow the patient's body weight to be absorbed proportionally
- 1.8 lb., 33 ILD upper layer and 1.8 lb., 40 ILD lower layer combine to provide the maximum pressure reduction and durability
- Reversible head-to-foot with contour cuts to reduce shear
- Anti-microbial waterproof cover with a 70 denier top to reduce friction and a vinyl bottom to protect the mattress
- Designed to meet EO184
- Does not contain latex



Invacare® CareGuard™ GM Gel Foam Mattress Overlay

Model no. CCGFMO

- Four large 2" thick gel bladders inside polyurethane foam
- Reversible for increased service life
- Durable aggressive base material with strong elastic corner securing straps reduce overlay movement
- Vapor permeable waterproof cover with low friction and shear (Do not machine wash)
- 250 lb. maximum patient weight
- Designed to meet EO185
- Does not contain latex



Invacare® CareGuard™ APP Alternating Pressure Pad System

Model no. CG9701

- Advanced pump technology alternately inflates and deflates air cells
- Variable pressure setting for added patient comfort
- Pump can be mounted on the bed with built-in bracket
- 250 lb. maximum patient weight
- Designed to meet EO180
- Does not contain latex

Specifications

Model no.	Weight	Length	Width	Thickness	Power	Code Approval	Warranty*	Flammability
CGSA1100	4 lb.	76"	33"	3"	n/a	n/a	6 months	n/a
CG10180	19 lb.	80"	35"	5"	n/a	n/a	2 years	Conforms to California Technical Bulletin #117
CCGFMO	55 lb.	76"	34"	3.5"	n/a	n/a	6 months	n/a
CG9701	Pump: 3.5 lb.	10"	6"	4"	115V, 160 HZ	UL-1431	2 years	
	Pad: 4.6 lb.	118"	34"	2.5"	n/a	n/a	30 days	n/a

* For manufacturer's defects