



**MODEL 600010 STRETCHER/CADAVER CARRIER**



Designed for multi-purpose use in hospitals, coroner's facilities and mortuaries. The Model No. 600010 is a durable and economical stainless steel stretcher/cadaver carrier.

**STANDARD DESIGN FEATURES**

- 1-1/4 inch tubular stainless steel construction.
- Heliarc welded for added strength and durability; ground and polished to a smooth finish.
- Four heavy duty 8 inch casters for easy movement (All four with brake mechanism)
- 18 gauge stainless steel top (With or without open corners)
- Top is bolted to lower frame to allow for removal, yet remains sturdy while in use
- Heavy-duty rubber bumper around perimeter.
- Dimensions 72 inches long x 26 inches wide x 34 inches high.