



MODEL 600039 COVERED CADAVER CARRIER



The COVERED CADAVER CARRIER allows for the concealed removal and transport of a deceased patient from the hospital room to the morgue. The carrier gives the illusion of being a standard empty gurney while hiding the cadaver beneath a specially fabricated top cover. The unit features stainless steel construction and adjustable elevation via a hand crank mechanism.

STANDARD FEATURES:

- Heavy gauge stainless steel frame and cadaver tray
- Sturdy, lightweight aluminum top frame to conceal cadaver - easily removable
- Custom fabricated cover made of black out material
- Hand crank mechanism elevates cadaver tray to 38-1/2 inches
- 8 inch casters, all with brake mechanism

Options: One location Total Wheel Lock
 Lever controlled directional caster for ease in steering