

RECONDITIONED BEDS

Reconditioned bed inventory includes Med/Surg, ICU and Birthing Beds.

All beds are visually inspected for bent or broken parts. Any bent or broken parts are repaired or replaced so that operating features function properly.

Our beds are completely disassembled. Our beds are cleaned inside and outside of the frame and all stickers or decals are removed. Our bed frames, arm rails and pans (if applicable) are sanded completely. All beds and parts are then cleaned again to remove any dust or dirt. Bed frames, arm rails and pans are painted in our professional paint booth. All parts and frames dry for 24 hours. Beds are then reassembled and include the following NEW parts:

- NEW OEM casters
- NEW plastic parts on arm rails
- NEW head and foot boards
- NEW labels or decals
- NEW capacitor and lubricant for motors

A final inspection of the beds operating features is performed and it is tested for electrical safety to ensure performance meets manufacturers' standards. Beds are blanket wrapped, covered in protective plastic and positioned in a custom stand for stability during shipment. Head and foot boards are shipped unattached and separately. Our current inventory includes: Hill-Rom Advanced, Versacare, Advanta, Total Care and Stryker Secure 3002 EX, Epic II, MPS beds.

Warranty: One year replacement parts warranty with lifetime technical assistance by phone.

