



MD101 Clinical Wall Hugger

Save Space Where You Need It Most



MD101 Clinical Wall Hugger

Standard Features:

- Mechanism moves forward as seat reclines. Seat can be placed close to wall, saving valuable floor space
- Hardened, Powder Coated, Steel Frame and Platform
- 1" Stainless Steel Push Handle
- Four H/D Casters
(Locking, Swivel Casters in Rear)
- Easily Removable/ Replaceable Components and Slip covers
- Outside Cable Release Mechanism
- High Density Foam is fire retardant and anti-microbial
- Vast variety of medical vinyls and fabrics available. Two-tone at no additional charge.



Standard Benefits:

- Exceptional Comfort for Patient
- Exceptional Durability and Cleanability for Caregiver
- Style enhances Clinical Décor

Optional Features:

- Four locking casters
- Fold Down Table (Left / Right)
- I V Pole with Brackets

Scope MD Seating Solutions for Quality Care