



# MD103 Clinical Recliner

## With Trendelenburg

### MD103 Clinical Recliner With Trendelenburg



#### Standard Features:

- True Trendelenburg
- Hardened, Powder Coated, Steel Frame and Platform
- 1" Stainless Steel Push Handle
- Four H/D Locking Casters (Swivel Casters in Rear)
- Swing Out Arms (Left & Right)
- Easily Removable/ Replaceable Components and Slip covers
- Outside Cable Release Mechanism
- High Density Foam is fire retardant and anti-microbial
- Vast variety of medical vinyls and fabrics available. Two-Tone at No Additional Charge

#### Standard Benefits:

- Exceptional Comfort for Patient
- Exceptional Durability and Cleanability for Caregiver
- Style enhances Clinical Décor

#### Optional Features:

- Fold Down Table (L, R, Both)
- I V Pole w/ Brackets

**Scope MD** ..... Seating Solutions for Quality Care