

HISTORY

Optical instruments have been a Seiler family tradition since 1913 when company founder Eric. H. Seiler entered the Zeiss-Jena School of Fine Optics in Jena, Germany. With a Master's degree in Fine Optics, Eric H. Seiler brought his craft to the United States where he and his wife, Dora L. Seiler, started Seiler Instrument and Manufacturing Company in 1945. His appreciation of fine craftsmanship and pride in his work contributed toward a philosophy aimed at quality and customer service, a legacy he passed on to his son, current president, Eric P. Seiler.

With Eric P. Seiler's business management skills and his father's extensive knowledge of optics, the company flourished and outgrew its original downtown St. Louis office. Seiler is now located in a suburban plant especially designed for optical instrument production. Today, a third generation has joined the company, sons Rick and Tom Seiler, whose management and sales skills contribute toward Seiler's growth and success.

Seiler keeps abreast of new technology with more sophisticated tools being added every year. Technology may change, but one thing continues to endure – Seiler's dedication to excellence and customer service, a philosophy that breeds customer satisfaction.

The use of Surgical Loupes offer many benefits. The obvious benefit is better vision. With magnifications up to 3.5X, the surgical loupes provide vision not obtainable with the naked eye. This improved vision allows you to perform most medical procedures more effectively. As the general public becomes more aware of this improved treatment, it is just a matter of time before your patients will demand the use of loupes.

FEATURES

- **2.5x, 3.0x and 3.5x** magnifications. Provides vision not attainable with the naked eye.
- **HIGH RESOLUTION** – highest grade glass for excellent clarity.
- **WATERSEALED** – optics are sealed to prevent dust and moisture.
- **DOUBLE HINGE FLIP-UP DESIGN** – for superior optical alignment.
- **LIGHTWEIGHT** – All frames are 100% titanium.
- **GREAT DEPTH OF FIELD** – comfortably extends the natural accommodation capability of the human eye.
- **WIDEFIELD** – highly corrected optics that allows the entire mouth to be viewed.
- **WORKING DISTANCE** – Short (340mm or 13.4 inches), Regular (420mm or 16.5 inches), Long (500mm or 19.7 inches).

All Seiler loupes allow the user to set your intrapupillary, convergence and declination angle for perfect optical alignment.

SERIES	250	300	350
Magnification (approx.)	2.5x	3.0x	3.5x
Working Distance			
Short	13.4"	13.4"	13.4"
Regular	16.5"	16.5"	16.5"
Long	19.7"	19.7"	19.7"

AVAILABLE COLORS

Metallic
Blue



Metallic
Red



Black



Metallic
Silver



2.5x

250S	2.5X 340mm	\$495. ⁰⁰
250R	2.5X 420mm	\$495. ⁰⁰
250L	2.5X 500mm	\$495. ⁰⁰



3.0x

300S	3.0X 340mm	\$545. ⁰⁰
300R	3.0X 420mm	\$545. ⁰⁰
300L	3.0X 500mm	\$545. ⁰⁰



3.5x

350S	3.5X 340mm	\$595. ⁰⁰
350R	3.5X 420mm	\$595. ⁰⁰
350L	3.5X 500mm	\$595. ⁰⁰

