



SURFACES

# Geo-Mattress with Wings<sup>®</sup>

**Category:** Therapeutic Raised Perimeter Mattress



## We Added Safety — Without Compromise!

The supportive, raised perimeter of our popular Geo-Mattress with Wings can help prevent resident injury. But instead of the same exposed, hard bolsters others use, its raised perimeter features a tissue-friendly Geo-Matt surface from edge to edge. Boasting all the key benefits of the *Max*, *Plus*, and *Pro*, Geo-Mattress with Wings means safety for your residents — and their skin!

### QUICK FACTS

#### Geo-Mattress<sup>®</sup> with Wings<sup>™</sup>

**Primary Use:**

- Minimize fall risk due to independent exits from bed

**Patient Indications:**

- General patient population, including:
- Moderate – High Risk for skin breakdown (e.g., Braden score  $\geq 10$ )
- Decreased mobility
- Repositions frequently by self or caregiver

**Key Advantages:**

- Edge-to-edge antishearing surface
- Extra-firm edge bolster
- Shearing protection at transfer opening
- Heel Slope<sup>®</sup>

#### Item # Description

W8435-13	Geo-Mattress with Wings, 84”L x 35”W x 6”H center, 8”H side
W8035-13	Geo-Mattress with Wings, 80”L x 35”W x 6”H center, 8”H side
W7535-13	Geo-Mattress with Wings, 75”L x 35”W x 6”H center, 8”H side

**Weight limit:** 400 lbs.  
**Warranty:** 5 years, non-prorated  
**Mattress Weight:** 24 lbs.

**Fire barrier versions** can be ordered to conform to NFPA 101 (Life Safety Code); Cal Tech Bulletins # 117, 129, and 603; 16 CFR Parts 1632 and 1633; and BFD IX-11. To order, replace "-13" suffix with "-29" suffix in model number.

# Geo-Mattress Atlas<sup>®</sup>

**Category:** Bariatric Therapeutic Foam Mattress



## Heavy Duty Protection for Your “Heftiest” Skin Integrity Challenges

The best of the clinically proven Geo-Mattress series, now available to meet the unique pressure management challenges of the overweight patient.

#### Item # Description

A7535-13	Geo-Mattress <sup>®</sup> Atlas, 75”L x 35”W x 7”H
A8035-13	Geo-Mattress <sup>®</sup> Atlas, 80”L x 35”W x 7”H
A8435-13	Geo-Mattress <sup>®</sup> Atlas, 84”L x 35”W x 7”H
* A8039-13	Geo-Mattress <sup>®</sup> Atlas, 80”L x 39”W x 7”H
* A8042-13	Geo-Mattress <sup>®</sup> Atlas, 80”L x 42”W x 7”H
* A8048-13	Geo-Mattress <sup>®</sup> Atlas, 80”L x 48”W x 7”H
* A8054-13	Geo-Mattress <sup>®</sup> Atlas, 80”L x 54”W x 7”H

**Additional bariatric options available – see page 15.**

**\*Allow up to 10 working days lead time for these sizes.**

**Weight limit:** 600 lbs. (35”); 675 lbs. (39”); 750 lbs. (42”, 48” and 54”)  
**Warranty:** 3 years, non-prorated  
**Mattress Weight:** 26 lbs. (35”)

**Fire barrier versions** can be ordered to conform to NFPA 101 (Life Safety Code); Cal Tech Bulletins # 117, 129, and 603; 16 CFR Parts 1632 and 1633; and BFD IX-11. To order, replace "-13" suffix with "-29" suffix in model number.

### QUICK FACTS

#### Geo-Mattress<sup>®</sup> Atlas

**Primary Use:**

- Prevention of skin breakdown

**Patient Indications:**

- Patient is 350-600 lbs.
- At risk for skin break down

**Key Advantages:**

- Progressive resistance accommodates wide range of patient weights
- Geo-Matt<sup>™</sup> anti-shearing surface
- Heel Slope<sup>®</sup>