

# PressureGuard Renew®

**Category:** Static Air/Foam Flotation

## Make Therapeutic Support a Renewable Resource

No matter how effective, typical all-foam mattress replacements share one common characteristic: their first day is their best day. While typical mattresses may continue to deliver acceptable pressure dispersion throughout their life, the orthopedic support they offer – especially at the lower back – inevitably degrades over time.

But the PressureGuard Renew is different.

Unlike all-foam and innerspring mattress replacements, fatigue won't rob it of its support and comfort over time.


Instead, the Renew maintains like-new performance both in pressure management and user support throughout its life. All it requires is a simple, occasional air system renewal that can be performed in seconds without ever removing the mattress from the bed.

Best of all, the Renew is equipped with all of the patented stability, safety, and simplicity innovations built into every PressureGuard model.



**Like-New Support in Seconds:**

*Add air through dual quick-connect hoses. Auto-adjust valve resets system to perfect level.*



### QUICK FACTS

#### PressureGuard Renew®

**Primary Use:**

- Prevention of skin breakdown

**Indications for Use:**

- Patient at moderate to high risk for skin breakdown (e.g. Braden score 12 or less)

**Key Advantages:**

- Proven effectiveness for skin management
- Refill air for like-new performance
- PressureGuard® design: stability, safety, anti-shearing

Item #	Description
2084-13	PressureGuard Renew®, 84"L x 35"W x 7"H
2080-13	PressureGuard Renew®, 80"L x 35"W x 7"H
* 3280-13	PressureGuard Renew®, 80"L x 32"W x 7"H
* 2078-13	PressureGuard Renew®, 78"L x 35"W x 7"H
2075-13	PressureGuard Renew®, 75"L x 35"W x 7"H
* RW7539-13	PressureGuard Renew®, 75"L x 39"W x 7"H <span style="float: right;">Twin Size</span>
* RW7554-13	PressureGuard Renew®, 75"L x 54"W x 7"H <span style="float: right;">Full Size</span>
* RW8060-13	PressureGuard Renew®, 80"L x 60"W x 7"H <span style="float: right;">Queen Size</span>
SERVICE-00	Air System Renewal Kit (includes hand pump and auto-adjust valve)

\* Allow up to 10 working days lead time for these sizes.

**Weight limit:** 500 lbs.  
**Warranty:** 5 years, non-prorated  
**Mattress Weight:** 20 lbs.

**Fire barrier versions** can be ordered to conform to NFPA 101 (Life Safety Code); Cal Tech Bulletins # 117, 129, and 603; 16 CFR Parts 1632 and 1633; and BFD IX-11. To order, replace "-13" suffix with "-29" suffix in model number.