

## RL-150 Rhinolaryngoscope\*

*Allows thorough exams of the upper respiratory system in less than five minutes.*

The RL-150 leverages Welch Allyn's vast experience in fiber optic imaging and distinguishes itself from other systems with the illumination of our patented Solarc™ Light Source. Delivers whiter, brighter light in a smaller, space-saving package.

### Features

- The thin (3.5 mm) insertion tube greatly reduces patient discomfort and will provide access to virtually all anatomical areas.
- Designed to complement your existing equipment, the RL-150 is adaptable to most light sources.
- The specifications include: 75° angle of view, 3-50 mm depth of view, 130°/130° deflection (up/down), 3.4 mm bending section diameter, 3.5 mm insertion tube diameter, and 30 cm insertion tube length. The RL-150 is immersible.
- For easy disinfection, the RL-150 Rhinolaryngoscope features an inexpensive disinfection tray that will reduce the amount of disinfectant used, thus minimizing the release of irritating fumes.
- The RL-150 provides a thorough exam of the upper respiratory tract in less than five minutes. The anatomy is easy to learn and the pathology is easily recognizable.
- Free *Primary Guide to Rhinolaryngoscopy* available.

### Model Numbers

65053: RL-150™ Rhinolaryngoscope and Solarc™ Light Source

65050: RL-150™ Rhinolaryngoscope Only

### Accessory Numbers

49500: Solarc™ Light Source, ACMI only

65002: RL-150™ Disinfection Tray

65003: RL-150™ Instrument Hanger

65004: Olympus CLK-3 Adapter

65005: Olympus CLK-4 Adapter

65009: Leakage Tester

65010: Olympus ILK-3 Adapter

65011: Machinda Adapter

65014: Pentax Adapter

65016: Wolf Adapter

65017: ACMI Adapter

65018: Storz Adapter

### Warranty

Welch Allyn will repair or replace, free of charge, any parts of its own manufacture proven to be defective through causes other than misuse, neglect, damage in shipment, or normal wear.


### Approvals

Light conforms to UL 2601-1; CAN/CSA C22.2 No. 601.1

\*Offered only in U.S.A. and selected countries.



RL-150 Rhinolaryngoscope, shown with 65002 Disinfection Tray



RL-150 Rhinolaryngoscope, shown with 65003 Storage Hanger