

18 767 Integrated Wall Systems

EYE, EAR, NOSE AND THROAT

767 Integrated Diagnostic Systems and Wall Transformer Sets

A convenient, modular package that puts everything needed for fast, accurate diagnoses at arm's length. By combining a wall aneroid, ophthalmoscope, otoscope, specula dispenser, and thermometer (optional), Welch Allyn puts all the basic diagnostic tools within easy reach. These systems are attractive, modular, and easy to order and install.



- > All the instruments you'll need to perform a complete diagnostic examination
- > Handle-based rheostats for quick light-intensity adjustments
- > Comfortable, contoured handles are easy to hold and maneuver
- > Durable, easy-to-clean, high-impact ABS housing
- > OptiSense™ sensor automatically turns handles on when they are removed from the cradles and off when they are replaced
- > Optional 200-specula KleenSpec® Plus Diagnostic or Pneumatic Specula Dispenser
- > Optional digital clock with pulse timer
- > Choose one of four halogen ophthalmoscopes, including the innovative PanOptic™ Ophthalmoscope that makes it easy to view the fundus 5X larger than ever before. Choose either PanOptic model 11810 or model 11820 with cobalt blue filter and corneal viewing lens. Other options include the traditional Standard Ophthalmoscope (11710) and the sealed (dust-free) Coaxial Ophthalmoscope (11720)
- > Otoscope options include the MacroView™, which provides 30% more magnification and nearly twice the field of view of traditional otoscopes and an ultrasealed Pneumatic Otoscope
- > 767 Wall Aneroid with easy-to-read dial face, certified accuracy and lifetime calibration program
- > Sturdy, predrilled mounting board mounts easily to any wall, preventing loss or theft
- > Optional 767 SureTemp® Electronic Thermometer for highly accurate oral readings in approximately 4-6 seconds, rectal readings in approximately 15 seconds

(See next page for full list of configurations)



76793-M
767 Integrated Diagnostic System with Coaxial Ophthalmoscope, MacroView Otoscope, Tycos Wall Aneroid, SureTemp Wall Thermometer and Otoscope Specula Dispenser

WELCH ALLYN INTEGRATED DIAGNOSTIC SYSTEM AND WALL TRANSFORMER SETS—COMPONENTS CHART












Part Number	MidMark Wall Board	767 Wall Transformer		Standard Ophthalmoscopes		PanOptic Ophthalmoscopes		Diagnostic MacroView Otoscopes		Pneumatic Otoscopes	767 Aneroid	KleenSpec Dispenser-Small	KleenSpec Dispenser-Large	SureTemp Thermometers	
	76790-1	76710	76720 (clock)	11710	11720	11810	11820	23810	23820	20200	7670-01	52401	52101	01767-200	01767-600 (w/ backlight)
76791-M	●	●		●				●			●	●		●	
76791-1MP	●	●				●		●			●	●		●	
76791-2MP	●	●					●	●			●	●		●	
76791-22MP	●		●				●	●			●	●		●	
76792-M	●	●		●				●			●	●			
76792-MNOBP	●	●		●				●			●	●			
76792-1MP	●	●				●		●			●	●			
76792-2MP	●	●					●	●			●	●			
76793-M	●	●			●			●			●	●		●	
76794-M	●	●			●			●			●	●			
76795-1MP	●	●				●			●		●	●		●	
76795-2MP	●	●					●		●		●	●		●	
76796-M	●	●			●				●		●	●			

20 767 Integrated Wall Systems

EYE, EAR, NOSE AND THROAT

WELCH ALLYN INTEGRATED DIAGNOSTIC SYSTEM AND WALL TRANSFORMER SETS—COMPONENTS CHART

Part Number	MidMark Wall Board	767 Wall Transformer		Standard Ophthalmoscopes		PanOptic Ophthalmoscopes		Diagnostic MacroView Otoscopes		Pneumatic Otoscopes	767 Aneroid	KleenSpec Dispenser-Small	KleenSpec Dispenser-Large	SureTemp Thermometers	
	76790-1	76710	76720 (clock)	11710	11720	11810	11820	23810	23820	20200	7670-01	52401	52101	01767-200	01767-600 (w/ backlight)
76796-1MP	•	•				•			•		•	•			
76796-2MP	•	•					•		•		•	•			
76797	•	•		•						•	•		•	•	
76797-2P	•	•					•			•	•		•	•	
76798-2P	•	•					•			•	•		•		
76710-71M		•		•				•							
76710-72M		•			•			•							
76710-81M		•				•		•							
76710-82M		•					•	•							
76710-L		•													•
76710-LA		•									•				•
76710-LAK		•									•	•			•
76710-LAP		•									•		•		•
76710-LK		•										•			•
76710-LP		•											•		•
76710-T		•												•	
76710-TA		•									•			•	
76710-TAK		•									•	•		•	
76710-TAP		•									•		•	•	
76710-TK		•										•		•	
76710-TP		•											•	•	
76720-L			•												•
76720-LA			•								•				•

WELCH ALLYN INTEGRATED DIAGNOSTIC SYSTEM AND WALL TRANSFORMER SETS—COMPONENTS CHART															
	MidMark Wall Board	767 Wall Transformer		Standard Ophthalmoscopes		PanOptic Ophthalmoscopes		Diagnostic MacroView Otoscopes		Pneumatic Otoscopes	767 Aneroid	KleenSpec Dispenser-Small	KleenSpec Dispenser-Large	SureTemp Thermometers	
															
Part Number	76790-1	76710	76720 (clock)	11710	11720	11810	11820	23810	23820	20200	7670-01	52401	52101	01767-200	01767-600 (w/ backlight)
76720-LAK			●								●	●			●
76720-LAP			●								●		●		●
76720-LK			●									●			●

Wall Set Accessories

Customize your Wall Set to match your individual needs.

76710	767 Wall Transformer	52101	KleenSpec® Pneumatic, Operating, and Consulting Oscope Disposable Specula Dispenser with Storage Compartment
76720	767 Wall Transformer with Digital Clock	52401	Universal KleenSpec® MacroView Oscope Disposable Specula Dispenser with Storage Compartment
76730	767 Third Handle Module	76700	Locking Collar for theft protection (set of 2)
76790-1	MidMark Wall Board		
01767-200	SureTemp® 4 Wall Thermometer		
01767-600	SureTemp® 4 Wall Thermometer with Backlight		
7670-01	767 Tycos® Wall Aneroid with Durable One-Piece Adult Cuff		