

Monitor any patient, any age, any location, any time.

Options and Accessories

The Propaq® Encore is a very intuitive monitor to use:

High resolution, bright EL display

This display, utilizing ICE technology, provides clear, bright, comprehensive patient data.

Quick-action hard keys

Fast operation of NIBP, alarms and ECG freeze functions to save you time!

Lightweight

Starting at 6.25 pounds, less than 13 pounds fully configured.

Rugged

Can withstand shock up to 50 g.

Longest battery life

Up to 5 hours depending on configuration and usage.

Small footprint

Needs little surface space, or just a bedrail.

Complete and ready to use


Comes with ECG and NIBP accessories and an AC power adapter.

The Propaq Encore is configurable with:

- Heart/pulse rate
- ECG (3 or 5 lead)
- Noninvasive blood pressure
- Invasive blood pressure (up to 2 channels)
- Temperature (2 channels)
- Nellcor® pulse oximetry
- Capnography (mainstream/sidestream)
- Impedance respiration
- Apnea



Propaq Encore allows you to bring the monitor to the patient, rather than move the patient to a fixed monitoring location.



View the most recent blood pressure (cuff) readings plus tabular trends and vital sign numerics.



Connect a Propaq Encore to the Acuity Central Monitoring Station as easily as you would plug in a phone.



The Modem Propaq uses standard telephone lines to transmit vital signs data to the Acuity Central Monitoring Station. From any location, patients can benefit from in-hospital surveillance and expert clinical analysis.




When you network your Propaq Encore monitors to the Acuity Central Monitoring Station, you'll see a map of your department, including up to 60 monitored patients, their status and location.

Linked to the Acuity® Central Monitoring Station, the Propaq Encore is the key to the Flexible Monitoring® system.

The Propaq Encore provides more monitoring for more patients for less money. Simply move it from bed to bed as your monitoring needs change. That's Flexible Monitoring, a proven clinical solution that reduces both acquisition and lifetime costs while better managing the process of care.

By listening to our customers, we have developed solutions that are saving our average hospital partner more than \$1,000,000 in operating costs per year. Just ask for more information and published studies documenting these patient care cost savings.



Shown at actual size.

Printer Option

The Propaq Encore printer provides high resolution, fully annotated real-time printouts that simultaneously document up to three waveforms and numerics for all active vital signs.

NIBP TREND 06/10/95						
TIME	HR/PR	SpO2	SYS / DIA	- MEAN	RR	
HH:MM	BPM	%	mmHg		RPM	
11:38	51	99	117 / 65	77	9	
11:37	60	99	119 / 66	81	14	
11:31	54	99	112 / 64	76	9	
11:30	56	99	114 / 64	76	8	
11:29	58	99	112 / 67	80	19	
11:29	58	99	112 / 65	79	12	
11:28	57	99	122 / 57	76	16	
11:26	59	98	122 / 66	78	17	

NIBP tabular trend printouts can provide up to 128 readings.



Mounting Systems

Tilt/swivel or fixed wall mounts, pole mounts, and roll stands, all with accessory baskets, can be added to optimize usage of your Propaq Encore monitors. Manufactured by GCX Corporation, these mounting systems are detailed in a separate brochure.

Carrying Cases

Heavy-duty Cordura® nylon carrying cases are available to protect your new Propaq Encore. There are five case styles, each one customized to a specific set of options and accessories.

