

#313/#342 Steel Frame Privacy Screens



ISO 9001-2000 CERTIFIED

Style, Function, and Versatility at an Affordable Price

Winco's steel frame privacy screens stand up to years of use. The wide, sturdy, steel frame eliminates tipping and leaning, while the twin wheel casters make maneuvering a breeze. Our

#313 screen can even be converted to a crutch tip base.



#342
Single Panel
Privacy Screen



#313
Three Panel
Folding Privacy Screen

#313/#342 Steel Frame Privacy Screens



ISO 9001-2000 CERTIFIED

Hinged Panel Ensures Years of Problem-Free Use

Flame retardant bacteriostatic vinyl curtains resist stains and are easily cleaned.



Screen Folds for Space-Saving or Storage

Smooth, Powder-Coated Steel Frame for Durability

Vinyl Curtains can Easily be Replaced by Removing the Spring Loaded Rods

Smooth Rolling Twin Wheel Casters (Includes Crutch Tips)



Wide Base for Added Stability



#313 SPECIFICATIONS

#342 SPECIFICATIONS

BASE SIZE:	20" x 27" (50.8 x 68.9cm)
OVERALL HEIGHT:	69" (175.3cm)
OVERALL WIDTH:	82" (208.3cm)
SINGLE PANEL WIDTH:	27" (68.6cm)
FLOOR TO BOTTOM OF CURTAIN:	17.5" (44.5cm)
FLOOR TO TOP OF LEG:	14.75" (37.5cm)
FRAME TYPE:	Powder Coated Steel
WEIGHT:	25lbs. (11.3kg)

BASE SIZE:	20" x 27" (50.8 x 68.9cm)
OVERALL HEIGHT:	69" (175.3cm)
OVERALL WIDTH:	42" (106.7cm)
SINGLE PANEL WIDTH:	42" (106.7cm)
FLOOR TO BOTTOM OF CURTAIN:	17.5" (44.5cm)
FLOOR TO TOP OF LEG:	14.75" (37.5cm)
FRAME TYPE:	Powder Coated Steel
WEIGHT:	15lbs. (6.8kg)