


#317 Designer Privacy Screen



ISO 9001-2000 CERTIFIED

Coordinates with the Other Winco Products in Your Facility

Winco's newest Privacy Screen has all the features you need, the reputation you trust, and a contemporary design. Match the curtain color to any standard Winco vinyl color, and you have a custom designed privacy screen that coordinates with the other Winco products in your facility!



#317
Three-Panel, Folding
Designer Privacy Screen

Proudly Made in the USA!



#317 Designer Privacy Screen



ISO 9001-2000 CERTIFIED

Hinged Panel Ensures Years of Problem-Free Use


Smooth, Powder-Coated Steel Frame for Durability

High Grade Anti-Bacterial, Anti-Microbial, Stain-Resistant Vinyl

Quick Assembly Legs - No Tools Required!

Smooth-Rolling Twin Wheel Casters

Wide Base for Superior Stability



VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors.)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Black	White	Burgundy	Mauve	Blueridge	Grotto
02	03	04	06	07	08	09	14	16	17	32
OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
33	34	37	41	42	43	44	60	61	62	63
Eucalyptus	Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress

BASE SIZE:	25" x 31" (63.5cm x 78.7cm)
OVERALL HEIGHT:	70" (177.8cm)
OVERALL WIDTH:	84" (213.4cm)
SINGLE PANEL WIDTH:	26" (66cm)
FLOOR TO BOTTOM OF CURTAIN:	13" (33 cm)
FLOOR TO TOP OF LEG:	8.5" (21.6cm)
FRAME TYPE:	Powder Coated Steel
WEIGHT:	45.5 lbs. (20.6 kg)