

# #653N/#654N Care Cliners



ISO 9001-2000 CERTIFIED

## Outstanding Value and Quality In Our Standard Care Cliner

Lower price does not mean lower quality when you buy a Winco! Our standard Care Cliner comes with all the superior features you might find on higher-priced clinical chairs.

It's no wonder Winco Care Cliners are among the bestselling clinical chairs on the market today!



*For Use In Dialysis,  
Oncology, Acute Care &  
Other Clinical Environments*

#653 N  
Standard Care Cliner  
w/ Nylon Casters



#654N  
XL Care Cliner  
w/ Nylon Casters



Proudly Made in the USA!

# #653N/#654N Care Cliners



ISO 9001-2000 CERTIFIED



Standard Trendelenberg Release on Back with Easy-Access Handle

Built-in, Ergonomic Headrest

Webbed Elastic and an Open Seat Base Provide Superior Comfort and Support

Smooth, Round Tube Frame With Recessed Hand-Grips

Adjustable Elastic Bands on Back Frame Provide Exceptional Comfort

Non-Porous, Fold-Down Side Tray Contains Small Spills - Standard on both sides

Removable Side-Panels for Easy Cleaning and Maintenance

High-Impact, Nylon, 5" (12.7cm) Casters  
Directional Lock on Left Rear  
All Others are Total Lock

## VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors.)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
<b>02</b>	<b>03</b>	<b>04</b>	<b>06</b>	<b>07</b>	<b>14</b>	<b>16</b>	<b>17</b>	OPTIONAL <b>32</b>	OPTIONAL <b>33</b>
OPTIONAL <b>34</b>	OPTIONAL <b>37</b>	OPTIONAL <b>41</b>	OPTIONAL <b>42</b>	OPTIONAL <b>43</b>	OPTIONAL <b>44</b>	OPTIONAL <b>60</b>	OPTIONAL <b>61</b>	OPTIONAL <b>62</b>	OPTIONAL <b>63</b>
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress







### #653N SPECIFICATIONS

### #654N SPECIFICATIONS

### OPTIONS

SEAT HEIGHT: 21" (53.3cm)  
 SEAT WIDTH: 21" (53.3cm)  
 SEAT DEPTH: 20" (50.8cm)  
 OVERALL HEIGHT: 49" (124.5cm)  
 OVERALL WIDTH: 31" (78.7cm)\*  
 BACK HEIGHT: 31.5" (80cm)  
 WEIGHT: 94 lbs. (42.6 kg)  
 WEIGHT CAPACITY: 275 lbs. (124.7 kg)  
 FOOTPRINT: 31" x 39" (78.7cm x 99.1cm)\*  
 FOOTPRINT (RECLINED): 31" x 76.5" (78.7cm x 194.3cm)\*

SEAT HEIGHT: 21" (53.3cm)  
 SEAT WIDTH: 25" (63.5cm)  
 SEAT DEPTH: 20" (50.8cm)  
 OVERALL HEIGHT: 49" (124.5cm)  
 OVERALL WIDTH: 35" (88.9 cm)  
 BACK HEIGHT: 31.5" (80cm)  
 WEIGHT: 111 lbs. (50.4 kg)  
 WEIGHT CAPACITY: 450 lbs. (204 kg)  
 FOOTPRINT: 35" x 44" (88.9cm x 111.8cm)  
 FOOTPRINT (RECLINED): 35" x 76.5" (88.9cm x 194.3cm)

-  **TB-133 Where Required**  
#653NT #654NT
-  **Steel Casters**  
#653 #654
-  **Heat & Massage**  
#653NH #654NH
-  **Heat ONLY**  
#653NW #654NW
-  **IV Pole**  
#IH
-  **Headrest Cover**  
#HC #HX

\* ADD 9" (22.9cm) TO WIDTH PER SIDE TABLE WHEN RAISED.