

#656 Premier Care Cliner



ISO 9001-2000 CERTIFIED

Style, Comfort, Value, and Great Features Make Our Premier Care Cliner Unlike any Other

Our Premier Care Cliner has a uniquely smooth, contemporary design to help your patients feel at ease, while maintaining the quality, durability, and great features you've come to expect from a Winco.

It's no wonder Winco Care Cliners are among the best selling clinical chairs on the market today!



#656
Premier
Care Cliner

For Use in Dialysis,
Oncology, Acute Care,
& Other Clinical Environments




Proudly Made in the USA!

#656 Premier Care Cliner



ISO 9001-2000 CERTIFIED

Adjustable Elastic Bands on Back Frame Provide Exceptional Comfort



Elastic Open Seat Base Provides Superior Comfort and Support

Non-Porous, Fold-Down Side Tray Contains Small Spills
Standard on Right Side
Optional on Left Side



True Trendelenberg Positioning Comes Standard

Low-Profile, Total-Lock, Steel Casters Provide Smooth, Easy Maneuvering

VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors.)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
02	03	04	06	07	14	16	17	32	33
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress
34	37	41	42	43	44	60	61	62	63

#656N SPECIFICATIONS

SEAT HEIGHT:	20" (50.8cm)
SEAT WIDTH:	21" (53.3cm)
SEAT DEPTH:	19.5" (49.5cm)
OVERALL HEIGHT:	46.5" (118.1cm)
OVERALL WIDTH:	30.5" (77.5cm)*
BACK HEIGHT:	30" (76.2cm)
WEIGHT:	93 lbs. (42.2 kg)
WEIGHT CAPACITY:	275 lbs. (124.7 kg)
FOOTPRINT:	30.5" x 43" (77.5cm x 109.2cm)*
FOOTPRINT (RECLINED):	30.5" x 76.5" (77.5cm x 194.3cm)*

* ADD 10" (25.4cm) TO WIDTH PER SIDE TABLE WHEN RAISED.

OPTIONS



TB-133
Where Required
#656T



2nd Side Table On Left
#656D



IV Pole
#IM



Headrest Cover
#HC